

2025 EDITION

INTERPRETIVE GUIDE

GI-MAP® – Unparalleled DNA-Based Stool Testing Our mission: to deliver innovative, accurate and clinically relevant diagnostic testing in a timely and cost-effective manner



GI Microbial Assay Plus



THANK YOU FOR CONSIDERING US!

"At Diagnostic Solutions Laboratory, we're not content with the range of clinical testing currently available to practitioners. We believe that every patient should achieve optimal health, and we're driven to give clinicians the tools to do so. Our mission, therefore, is to use our resources to bring the most advanced, innovative, and clinically relevant testing to healthcare providers worldwide."

Tony Hoffman
President and CEO



CONTENTS

THANK YOU	2
INTRODUCTION	4
HOW TO READ THE REPORT	5
PATHOGENS	6
H. PYLORI AND VIRULENCE FACTORS	16
COMMENSAL/KEYSTONE BACTERIA	20
OPPORTUNISTIC BACTERIA	23
FUNGI/YEAST	27
VIRUSES	29
PARASITES	31
INTESTINAL HEALTH MARKERS	36
ADD-ON TESTS	40
ANTIBIOTIC RESISTANCE GENES	42
STOOLOMX [™]	44
GI-MAP PATTERNS	56
SUPPORT INFO	63



INTRODUCTION

The Gastrointestinal Microbial Assay Plus (GI-MAP®) is an innovative clinical tool that measures gastrointestinal microbiota DNA from a single stool sample with state of the art, quantitative polymerase chain reaction (qPCR or real-time PCR) technology.

The GI-MAP was designed to detect microbes that may be disturbing normal microbial balance or contributing to illness as well as indicators of digestion, absorption, inflammation, and immune function. The following guide may be useful for understanding the nature of each of the microorganisms found on the GI-MAP, as well as clinical implications and treatment guidelines.

Please see the GI-MAP white paper for additional, fully referenced, information.





HOW TO READ THE REPORT

GI-MAP quantifies bacteria, fungi, viruses, and parasites using qPCR. This is a leap forward from older methodologies that report only positive or negative. Results are reported as colony forming units per gram of stool (CFU/g). One CFU is roughly equivalent to one microorganism (or one cell). Results are expressed in standard scientific notation. A reported result of 3.5e7 is equivalent to 3.5×10^7 CFU/g, which equals 35,000,000 CFU/g, or 35 million CFU per gram of stool.

	PATHOG	ENS	
BACTERIAL PATHOGENS	Result		Reference
Campylobacter	< dl		< 1.00 e3
C. difficile Toxin A	1.21 e5	High †	< 1.00 e3

Figure 1. The normal reference range for *C. difficile*, Toxin A is 0–1,000 CFU/g. The patient's result is very high at 1.21×10^5 , or 121,000 CFU/g.

Reference ranges were developed using known positive, diseased samples to construct cut off values that distinguish disease-causing amounts of pathogenic and opportunistic microbes. Reference ranges for the pathogens were correlated with an FDA cleared assay for GI pathogens. The GI-MAP is capable of detecting as low as 0.1 cell per gram of stool.

Table 1. Scientific notation; a basic reference table.

1.0e1	1 X 10 ¹	10	Ten
1.0e2	1 X 10 ²	100	One hundred
1.0e3	1 X 10 ³	1,000	One thousand
1.0e4	1 X 10 ⁴	10,000	Ten thousand
1.0e5	1 X 10 ⁵	100,000	One hundred thousand
1.0e6	1 X 10 ⁶	1,000,000	One million



PATHOGENS

The GI-MAP includes pathogens (bacterial, parasitic, and viral) commonly known to cause intestinal gastroenteritis. It's important to note that not all individuals with positive findings for pathogens will present with symptoms. Many factors, including the health of the individual, the transient nature of some pathogens, and the presence and expression of virulence factors all contribute to an individual's symptoms. Toxins are a type of virulence factor produced by certain pathogens. Since GI-MAP is a DNA-based test, results reflect the levels of pathogenic strains carrying the toxin genes, not the levels of any toxins that may be produced.

BACTERIAL PATHOGENS

Campylobacter

- Epidemiology
 - » One of the most common causes of foodborne illness in the U.S.
 - » Fecal contamination of poultry and water
- Clinical Implications
 - » May be infectious at very low exposures
 - » Symptoms range from mild to severe abdominal pain, diarrhea,

- fever, malaise; lasting several days to several weeks
- » Vast majority of those with symptoms of gastroenteritis recover without treatment
- Therapeutic Approaches
 & Considerations
 - » See patient's calprotectin level to determine GI inflammation
 - » Consider high dose probiotics, broad-spectrum antimicrobial herbs, and 5R Protocol (see Table 2)
 - » Heavy infections can be treated with azithromycin and fluoroquinolones





Clostridium difficile, Toxin A and Toxin B

The GI-MAP tests only for the genes for toxin A and toxin B, which are carried by *C. difficile*. The GI-MAP does not measure toxins directly for any microbe.

Epidemiology

- » 2–10% of population are carriers, most are asymptomatic
- » Prolonged use of antibiotics may be causative factor

• Clinical Implications

- » Symptoms include inflammation, abdominal pain, cramping, fever, and diarrhea
- » Symptoms often present during antibiotic use and often subside once antibiotics are discontinued
- » Gastrointestinal (GI) infection can cause reactive arthritis

Therapeutic Approaches & Considerations

- » See patient's calprotectin and secretory IgA (SIgA) levels to determine GI inflammation and immune response.
- » Consider Saccharomyces boulardii, high dose probiotics, broadspectrum antimicrobial herbs, and 5R Protocol (see Table 2)
- » Mild infections can be treated with metronidazole
- » Heavy infections can be treated with vancomycin and fidaxomicin
- » Asymptomatic patients may not need treatment

- » In asymptomatic patients with positive toxins A and/or B, the genes are likely not "turned on," and thus not causing disease. It is still prudent to avoid antibiotics in these patients to prevent CDAD. Consider antimicrobial herbal formulas, which can suppress *C. diff* without activating toxin production.
- » Additional testing for toxins A and B may be warranted

E. coli EPEC/EHEC

Epidemiology

- » Fecal contamination of food (undercooked beef, raw milk, and unpasteurized juice) and water
- » Implicated in hemorrhagic colitis, may progress to hemolytic uremic syndrome (HUS)

Clinical Implications

- » Symptoms include fever, abdominal cramping, fatigue, nausea, and diarrhea
- » Symptoms may last up to a week

- » See patient's calprotectin and SIgA levels to determine GI inflammation and immune response.
- Antibiotics may be contraindicated;
 they can initiate hemolytic
 uremic syndrome (HUS)
- » Consider high-dose probiotics (300+ billion CFU/g) such as: *Lactobacillus acidophilus, Bifidobacterium bifidum, Bifidobacterium longum, Lactobacillus rhamnosus,*



Table 2. Clinical Approach — *The Five "R" Treatment Protocol.*

The 5R Protocol is a widely accepted clinical guideline to treating pathogens and imbalances in the GI microbiota and restoring health to the gastrointestinal tract. Re-test patients with the GI-MAP in 3–6 months to monitor progress and make changes to the treatment protocol as needed.

	Antimicrobial	Broad-spectrum antimicrobial herbs including: berberine, caprylic acid, garlic oil, oil of oregano, uva ursi, olive leaf extract
REMOVE Using a course of antimicrobial,	Antibiotics	Research the recommended antibiotic for the specific microbe present. Avoid medications to which the microbe is thought to have resistance.
antiviral, antifungal, or antiparasitic therapies in cases where these organisms are	Antifungal	Caprylic acid, garlic oil, oil of oregano, olive leaf extract
present. It may also be necessary to remove offending foods, gluten, or medication that may be acting as antagonists.	Antiparasitic	Black walnut, garlic oil, oil of oregano, Artemisia (wormwood), berberine, goldenseal, gentian root extract, quassia bark extract, citrus seed extract
	Antiviral	Olive leaf extract, purified silver, cat's claw, monolaurin, osha root (<i>Ligusticum porteri</i>), vitamin A, vitamin C, vitamin D, reishi mushrooms, <i>Echinacea</i> , zinc
REPLACE In cases of maldigestion or malabsorption, it may be necessary to restore proper digestion by supplementing with digestive enzymes.	Digestive support	Betaine hydrochloride, apple cider vinegar, herbal bitters, ox bile, lactase, pancreatic enzymes (amylase, lipase, protease), pepsin
REINOCULATE Recolonization with healthy, beneficial bacteria. Supplementation with probiotics,	Probiotics	Lactobacillus acidophilus, Bifidobacterium bifidum, Bifidobacterium longum, Lactobacillus rhamnosus, Bifidobacterium breve, Lactobacillus casei, Saccharomyces boulardii
along with the use of prebiotics helps re-establish the proper microbial balance.	Prebiotics	Beta-glucan, fiber, inulin, pectin, xylooligosaccharides, galactooligosaccharides, larch arabinogalactans
REPAIR	Immune Support	Colostrum, immunoglobulins, <i>S. boulardii</i>
Restore the integrity of the gut mucosa by giving support to healthy mucosal cells, as well as immune support.	Intestinal Barrier Repair	L-Glutamine, aloe vera extract, deglycyrrhizinated licorice, marshmallow root, okra, N-acetyl glucosamine, quercetin, <i>S. boulardii</i> , slippery elm, zinc carnosine, vitamin A, essential fatty acids, B vitamins
REBALANCE Address whole body health and lifestyle factors so as to prevent future GI dysfunction.	Support Consideration	Sleep, diet, exercise, and stress management





- Bifidobacterium breve, Lactobacillus casei, Streptococcus thermophilus
- » Consider bacteriophages, broadspectrum antimicrobial herbs, and 5R Protocol (see Table 2)

E. coli 0157

- Epidemiology
 - » Fecal contamination of food and liquids (dairy, undercooked beef, vegetables, juices)
 - » Implicated in many outbreaks and cases of bloody diarrhea and HUS
 - » High prevalence worldwide

Clinical Implications

- » Symptoms may include severe abdominal cramps and diarrhea
- Shiga toxins inhibit protein synthesis& elicit strong inflammatory response
- Therapeutic Approaches & Considerations
 - » See patient's calprotectin and SIgA levels to determine GI inflammation and immune response.
 - » Antibiotics may be contraindicated; they can initiate HUS
 - Consider high-dose probiotics (300+ billion CFU/d)
 - » Consider bacteriophages, broadspectrum antimicrobial herbs, and 5R Protocol (see Table 2)

SOURCES OF EXPOSURE AND RE-INFECTION

To effectively treat infections and prevent reinfection, exposure should be identified and eliminated. Most exposure to pathogens occurs via fecal-oral transmission, most often due to use of contaminated water sources or improper hand hygiene. This may include drinking contaminated water, eating raw foods washed in contaminated water or harvested (e.g. shellfish) in contaminated water, or improper handwashing.

To remove microorganisms from food, the FDA recommends first washing your hands, running cool water over fruits and vegetables, while rubbing or scrubbing, and then letting them dry out before eating. During treatment, consider all possible sources of fecal transmission: romantic partners, children (especially if in diapers or not toilet-trained), sheets, towels, water source to the home, etc...



Enteroinvasive E. coli / Shigella

- Epidemiology
 - » Fecal contamination of ingested foods
- · Clinical Implications
 - » Symptoms include diarrhea (with blood and/or mucus), vomiting, fever, chills, fatigue, and abdominal cramping
 - » Symptoms are generally self-limiting
 - » Gastrointestinal (GI) infection can cause reactive arthritis
- Therapeutic Approaches & Considerations
 - » See patient's calprotectin and SIgA levels to determine GI inflammation and immune response.
 - » Antibiotics may be contraindicated; they can initiate HUS
 - » Consider high-dose probiotics (300+ billion CFU/d)
 - » Consider bacteriophages, broadspectrum antimicrobial herbs, and 5R Protocol (see Table 2)

Enterotoxigenic E. coli LT/ST

- Epidemiology
 - » Fecal contamination of ingested foods (hamburger meat, unpasteurized milk, and contaminated water)
- Clinical Implications
 - » Symptoms include watery, bloody diarrhea

Therapeutic Approaches & Considerations

- » See patient's calprotectin and SIgA levels to determine GI inflammation and immune response
- » Antibiotics may be contraindicated; they can initiate HUS
- » Consider high-dose probiotics (300+ billion CFU/d)
- » Consider bacteriophages, broadspectrum antimicrobial herbs, and 5R Protocol (see Table 2)

Shiga-like Toxin E. coli stx1

- Epidemiology
 - » Most common cause of traveler's diarrhea
- Clinical Implications
 - » Diarrhea is the most common symptom
- Therapeutic Approaches
 & Considerations
 - » See patient's calprotectin and SIgA levels to determine GI inflammation and immune response
 - » Antibiotics may be contraindicated; they can initiate HUS
 - » Consider high-dose probiotics (300+ billion CFU/d)
 - » Consider bacteriophages, broadspectrum antimicrobial herbs, and 5R Protocol (see Table 2)





Shiga-like Toxin E. coli stx2

- Epidemiology
 - » Fecal contamination of ingested foods (undercooked meat, unpasteurized milk, juice, and water)
- Clinical Implications
 - » Symptoms may include severe abdominal cramps and diarrhea
 - » Toxins may elicit strong inflammatory response
- Therapeutic Approaches & Considerations
 - » See patient's calprotectin and SIgA levels to determine GI inflammation and immune response
 - » Antibiotics may be contraindicated; they can initiate HUS
 - Consider high-dose probiotics (300+ billion CFU/d)
 - » Consider bacteriophages, broadspectrum antimicrobial herbs, and 5R Protocol (see Table 2)
 - » Antibiotics and antidiarrheal medicines are contraindicated; they may increase the risk of developing HUS

Salmonella

- Epidemiology
 - » Fecal contamination of ingested foods (eggs, poultry, meat, unpasteurized milk, raw fruits, and vegetables)
 - » Exposure to pets (reptiles, amphibians, baby chicks)
- Clinical Implications
 - » May be asymptomatic

- » Symptoms include fever, vomiting, and severe diarrhea
- » Typically self limiting within seven days
- » GI infection can cause reactive arthritis and may be involved in ankylosing spondylitis
- » Systemic infections may require treatment with antibiotics
- Therapeutic Approaches & Considerations
 - » See patient's calprotectin and SIgA levels to determine GI inflammation and immune response
 - » Remove sources of infection
 - Consider high-dose probiotics (300+ billion CFU/d)
 - Consider broad-spectrum antimicrobial herbs and 5R
 Protocol (see Table 2) Antibiotics are contraindicated; they may cause relapse of infection

Table 3.

FOOD SOURCES OF SALMONELLA

Poultry and poulty products

Meat

Dairy

Raw, fresh, ready-to-eat produce such as:

- Tomatoes
- · Leafy greens
- Sprouts
- Berries
- Melons



Vibrio cholerae

Epidemiology

» Fecal contamination of ingested foods (raw shellfish) and often picked up during international travel

Clinical Implications

- » May be asymptomatic or cause mild symptoms
- » Severe infections present with profuse watery diarrhea ("rice-water stools"), vomiting, rapid heart rate, loss of skin elasticity, thirst, dry mucous membranes, low blood pressure, restlessness, or irritability

Therapeutic Approaches & Considerations

- » See patient's calprotectin and SIgA levels to determine GI inflammation and immune response.
- » Rehydration therapy
- » Zinc, especially in children
- » Consider probiotics, broadspectrum antimicrobial herbs and 5R Protocol (see Table 2)
- » Heavy infections may be treated with doxycycline

Yersinia enterocolitica

Epidemiology

» Fecal contamination of ingested foods and liquids (water, undercooked pork, meats, and dairy products)

Clinical Implications

- » Symptoms usually develop four to seven days after exposure and are self-limiting
- » Symptoms include water or bloody diarrhea, fever, vomiting, and abdominal pain (may resemble appendicitis)
- » Symptoms may mimic Crohn's disease
- May trigger autoimmune thyroiditis or inflammatory arthritis in susceptible individuals

- » Consider probiotics, broadspectrum antimicrobial herbs and 5R Protocol (see Table 2)
- » Heavy infections can be treated with doxycycline in combination with an aminoglycoside
- » Trimethoprim-sulfamethoxasole, chloramphenicol, and rifaximin may also be useful treatments





PARASITIC PATHOGENS

Cryptosporidium

Epidemiology

- » Fecal contamination of ingested foods and liquids (contaminated water and swimming pools, undercooked meat, and raw milk)
- » Common cause of traveler's diarrhea

Clinical Implications

- » Symptoms typically last 2–3 weeks and are self-limiting
- » If symptoms persist, look for sources of contamination, such as drinking water
- » Can cause reactive arthritis

Therapeutic Approaches & Considerations

- » May not require treatment
- » See patient's calprotectin and SIgA levels to determine GI inflammation and immune response
- » If necessary, consider anti-parasitic herbal treatments containing ingredients such as black walnut, garlic oil, oil of oregano, Artemisia (wormwood), berberine, goldenseal, gentian root extract, quassia bark extract, citrus seed extract
- Consider probiotics and 5R Protocol (see Table 2)
- » Search for and remove sources of fecal contamination
- » Heavy infections can be treated with nitazoxanide*

Entamoeba histolytica

Epidemiology

- » Fecal contamination of ingested foods or water
- » Pets may be a source of exposure
- » Sexual contact may be a source of exposure

Clinical Implications

- » Symptoms include diarrhea, fulminating colitis (resembling ulcerative colitis), and dysentery
- » Extreme cases may invade liver and lung tissues

- » See patient's calprotectin and SIgA levels to determine GI inflammation and immune response
- » Treatment may be indicated, even in asymptomatic carriers
- » Mild infections can be treated with lodoquinol, paromomycin, or diloxanide furoate*
- » Moderate to heavy infections can be treated with metronidazole or tinidazole, followed by iodoquinol or paromomycin*
- » If appropriate, consider anti-parasitic herbal treatments (see Table 2)
- » Consider probiotics and 5R Protocol (see Table 2)
- » Avoid reinfection by fecal contamination



Giardia

Epidemiology


- » Most commonly isolated protozoan worldwide
- » Found in outside water sources (lakes, streams, ponds) and can get past filtration systems
- » Carried by animals
- » Common in daycare workers

Clinical Conditions

- » May be asymptomatic, especially in patients with adequate levels of normal bacteria and SIgA
- » Symptoms include acute diarrhea, bloating, cramps, weight loss, intestinal malabsorption, and steatorrhea
- » Can cause urticaria or neurologic symptoms such as irritability, sleep disorder, or depression

- » May cause malnutrition and vitamin B12 deficiency
- » Can cause reactive arthritis

- » See patient's calprotectin and SIgA levels to determine GI inflammation and immune response
- » Infections can be treated with tinidazole, nitazoxanide, metronidazole, paramomycin, furazolidone, or quinacrine*
- » Consider probiotics and 5R Protocol (see Table 2) to repair and rebuild the gut mucosa
- * Additional information (dosing, efficacy, etc.) on pharmaceutical treatment for parasites may be found at www.cdc.gov/parasites/index.html and in the Physician's Desk Reference.



Giardia intestinalis protozoan





VIRAL PATHOGENS

Adenovirus 40/41

Epidemiology

- » Common cause of diarrhea in infants and children but can also affect adults
- » Mainly transmitted by fecal contamination (fecal-oral route)

Clinical Implications

- » Causative agents of gastrointestinal disease and gastroenteritis
- » Symptoms include fever and watery diarrhea, usually limited to 1–2 weeks
- » May also be present in the stool of asymptomatic carriers and may not require treatment

Treatment

- » Handwashing
- » Hydration
- » Antiviral herbs such as cat's claw, osha root, reishi mushrooms, vitamins A, C, and D, zinc, Echinacea
- » Address other imbalances on the GI-MAP and use 5R Protocol (see Table 2) to rebuild gut health and gut immunity

Norovirus GI/GII

Epidemiology

- » Fecal contamination of ingested foods and water
- » Common cause of stomach flu on cruise ships
- » Common cause of nonbacterial gastroenteritis and outbreaks in the world

Clinical Implications

- » Symptoms include nausea and vomiting, diarrhea, abdominal cramps, low-grade fever, muscle aches, fatigue, and headache
- » Generally short-lived, lasting about 24–72 hours

Treatment

- » Antivirals are not recommended
- » Supportive care for the gastric mucosa, hydration, and immuneboosting agents may be warranted
 - Handwashing
 - Hydration
 - Antiviral herbs such as cat's claw, osha root, reishi mushrooms, vitamins A, C, and D, zinc, Echinacea
 - Address other imbalances on the GI-MAP and use 5R Protocol (see Table 2) to rebuild gut health



H. pylori AND VIRULENCE FACTORS

Helicobacter pylori (H. pylori)

Recent studies have shown that nearly 50% of the world's population may harbor *H. pylori*. And, although many carriers are asymptomatic, *H. pylori* is known to have a causative role in ulcers, chronic gastritis, and stomach cancer. Additionally, in early phases of colonization, patients may experience hypochlorhydria followed by a change to hyper aciduria. Over time, additional *H. pylori* strains may colonize, including those with Virulence Factors and increased disease potential.

Epidemiology

» Fecal contamination, oral to oral, and family inter-infection are common modes of transmission

• Clinical Implications

- » Dyspepsia, abdominal pain, nausea, vomiting and chronic gastrointestinal symptoms
- » Peptic ulcers
- » May induce mucosal atrophy and metaplastic changes

Virulence Factors

Of the 50% of the population believed to be infected with *H. pylori*, only 2% develop gastric cancer. Positive *H. pylori* virulence factors on the GI-MAP represent the *genetic potential* for an *H. pylori* strain to cause pathology. For example, some clinicians may choose an aggressive treatment protocol for a patient with dyspepsia and a family history of gastrointestinal cancer, who shows elevated *H. pylori* and positive virulence factors. (*See Table 4*)

Treatment

- » Asymptomatic patients may not require treatment
- » Consider herbal formulas to eradicate or suppress *H. pylori*. Ingredients may include: deglycyrrhizinated licorice, mastic gum, methylmethionine sulfonium chloride, vitamin C, zinc carnosine, bismuth citrate, berberine, goldenseal, oil of oregano, grape extract, Chinese goldthread extract, yerba mansa extract
- » See pancreatic elastase-1 to determine if maldigestion and/or hypochlorhydria might be present
- » Consider high-dose probiotics and 5R Protocol (see Table 2)






- » Rebuild healthy gastric mucosa by reducing stress and giving soothing and healing agents such as glutamine, aloe, DGL, and vitamin A
- » Address dental hygiene; the mouth is a reservoir for H. pylori
- Consider sources of exposure, especially romantic partners or family members
- » Address other imbalances on the GI-MAP
- » For peptic ulcer disease, the firstline triple therapy (prescription) treatment includes a proton pump inhibitor, clarithromycin, and amoxicillin or metronidazole
- » Fluoroquinolones and tetracycline are used in second-line regimens against *H. pylori*
- » See the GI-MAP Antibiotic Resistance Genes results for *Helicobacter* when designing an antibiotic protocol

H. PYLORI ANTIBIOTIC RESISTANCE GENES

	Result	Reference
Clarithromycin	Positive	Negative
Genes associated with clarithron	nycin resistance	
A2142C	Absent	
A2142G	Absent	
A2143G	Present	

Figure 2. An example antibiotic resistance gene measured on the GI-MAP for *H. pylori*. Number and letter combinations are single nucleotide polymorphisms (SNPs), or gene targets, that are involved in *H. pylori* drug resistance. If any SNP is detected (*present*), then the *H. pylori* strain/s are resistant to that class of antibiotics. In this example, the patient's *H. pylori* is resistant to the clarithromycin class of antibiotics and it would be prudent to use a different antibiotic when tailoring a treatment protocol.

Note: In supporting individuals with H. pylori, consider the patient history, clinical symptoms, virulence genes, and the amount of H.pylori DNA detected in order to design a customized treatment plan.



Helicobacter pylori bacterium



Table 4. *H. pylori* virulence factors and disease associations.



Helicobacter pylori Virulence Factors

Quick Reference Guide

The virulence factor genes on GI-MAP® are found exclusively on the genome of *H. pylori*.

These genes code for proteins that will predispose one to more serious *H. pylori* infections. The chart below provides details of each virulence factors tested on the GI-MAP.

Gene Acronym	Gene Name	Genetic Characteristics	Associations with Disease
BabA	Blood Group Antigen Binding Adhesion	 Promotes DNA breakage in host cell Improves H. pylori adherence ("stickiness") to epithelial cells May promote other virulence factors, especially CagA 	May promote carcinogenesis
CagA	Cytotoxin Associated Gene A	 Promotes <i>H. pylori</i> adhesion and colonization Affects barrier function of gastric epithelial tight junctions Promotes loss of cell polarity Antagonizes VacA Evades the immune system and affects the activity of dendritic cells and B-cells. Considered part of the "pathogenicity island" which includes VirB and VirD virulence factors. This is a closely-associated group of genes that work synergistically and often transfer as a unit. 	Promotes carcinogenesis, strong association. Also associated with peptic ulcer disease
DupA	Duodenal Ulcer- Promoting Gene A	Promotes inflammation	Associated with duodenal ulcers, specifically
IceA	Induced by Contact with Epithelium A	 Transcription of this gene is only initiated after adhesion to the gastric epithelium Promotes inflammation and associated with elevated IL-8 	Associated with dyspepsia and gastric & duodenal ulcers NOT associated with gastric cancel
OipA	Outer Inflammatory Protein A	Promotes inflammationDrives IL-8 production	Associated with carcinogenesis and peptic ulcer disease
VacA	Vacuolating Toxin A	 Enters the host cell by endocytosis Affects mitochondrial function Disrupts tight junctions Causes a programmed necrosis by inducing the production of large vacuoles inside the host cells; inducing cellular swelling; disrupting cell barrier thus causing nutrient leakage Facilitates nutrient acquisition (iron, minerals, amino acids, etc.) Inhibits antigen presentation in vitro Antagonizes CagA 	Associated with gastric inflammation, peptic ulcer disease and gastric cancers
VirB & VirD		 Part of the CagA "pathogenicity island" Both genes can potentiate CagA virulence factor by aiding in its transmission to host epithelial cels In the absence of CagA, these virulence factors are unlikely to change clinical outcome of <i>H. pylori</i> infections. 	Evaluate next to CagA virulence factors. VirB & VirD, if positive, ca potentiate CagA virulence and clinical associations





- How To Target Treatments in the Presence of Virulence Factors:
 - » Generally, when virulence factors are present, the treatment goal will be to fully eradicate the *H. pylori* population. This can be confirmed by retesting the full GI-MAP, the pathogen panel, or the *H. Pylori* panel 4–6 weeks after completing treatment. The goal is to achieve a result of <dl on the retest. The exception to this may be</p>
- VirB and VirD if they are found in isolation (without CagA present).
- » Below is a chart of treatment considerations for each of the virulence factors. These would be used in addition to the standard treatments for *H. pylori* alone.

Table 4a. *H. pylori* virulence factors and treatment considerations.

Virulence Factor	Special Treatment Considerations
BabA	More aggressive treatment may be warranted; consider the use of adhesion
	inhibitions, particularly cranberry
CagA	Target inflammatory support, promote T-cell activity, consider curcumin,
	resveratrol/red wine, ginger, Nigella sativa, low salt diet
DupA	Consider the use of demulcents for mucosal protection
IceA	Inflammatory support, consider the use of adhesion inhibitors
OipA	Inflammatory support
VacA	Mitochondrial support, consider Nigella sativa, green tea, red wine/resveratrol,
	Scutellaria baicalensis
VirB & VirD	No additional treatments necessary


Please see the GI-MAP white paper for additional, fully referenced, information.





COMMENSAL/ KEYSTONE BACTERIA

Trillions of microorganisms inhabit the human intestine to make up a complex ecosystem that plays an important role in human health. Commensal bacteria extract nutrients and energy from our diets, maintain gut barrier function, produce vitamins (biotin and vitamin K), and protect against colonization by potential pathogens.



Bacteria Lactobacillus, lactic acid bacteria which are part of normal flora of human intestine





THE FOLLOWING COMMENSAL/KEYSTONE BACTERIA ARE REPORTED ON THE GI-MAP

Table 5.

COMMENSAL BACTERIA	
Bacteroides fragilis	Gram-negative species of the <i>Bacteroidetes</i> phylum. Immune-modulating normal gut species. Believed to be involved in microbial balance, barrier integrity, and neuroimmune health (Hsiao 2013). High levels may result from reduced digestive capacity or constipation. Low levels may contribute to reduced anti-inflammatory activity in the intestine.
Bifidobacterium spp.	Gram-positive genus in the <i>Actinobacteria</i> phylum. Present in breast milk. Colonizes the human GI tract at birth. Common in probiotics. Thrives on a wide variety of prebiotic fibers. Low levels may result from low fiber intake or reduced mucosal health. High levels are more common in children than in adults.
Enterococcus spp.	Gram-positive genus of lactate-producing bacteria in the <i>Firmicutes</i> phylum. High levels may be due to reduced digestive capacity, constipation or small intestinal bacterial overgrowth. Low levels may indicate insufficiency of beneficial bacteria.
Escherichia spp.	Gram-negative genus in the <i>Proteobacteria</i> phylum. Normal gut flora. <i>Escherichia coli</i> (<i>E. coli</i>) is the primary species in this genus. Most <i>E. coli</i> are nonpathogenic (pathogenic <i>E. coli strains are measured separately in "Pathogens" section of the GI-MAP</i>). High levels may be indicative of increased intestinal inflammatory activity. Low levels may indicate reduced mucosal health and decreased protection against pathogenic <i>E. coli</i> .
Lactobacillus spp.	Gram-positive genus of lactate-producing bacteria in the <i>Firmicutes</i> phylum. Many strains used as probiotics. High levels may result from reduced digestive capacity or excessive intake of carbohydrates. Low levels may be due to low carbohydrate intake or high salt intake, and may also indicate reduced mucosal health.
Enterobacter spp.	Gram-negative genus in the <i>Proteobacteria</i> phylum. Closely related to <i>E. coli (in the same taxonomic family)</i> . High levels may indicate increased intestinal inflammatory activity. Low levels may indicate reduced mucosal health.
Akkermansia municiphila	Keystone species and primary mucus degrader. Generates mucus-derived sugars and metabolic products that support the growth and energy needs of other gut microbes. Promotes mucosal health and mucus production. Low levels associated with obesity and metabolic dysfunction. High levels linked to multiple sclerosis.
Faecalibacterium prausnitzii	Widely recognized as an important keystone species in the Clostridia class, as well as a major butyrate producer. Promotes anti-inflammatory processes and mucosal homeostasis. Reduced levels have been associated with a wide range of chronic inflammatory and autoimmune diseases.
Roseburia spp.	A genus of Gram-positive anaerobic bacteria in the <i>Clostridia</i> class that inhabit the human colon. The <i>Roseburia</i> genus has five well-characterized species, all of which produce short chain fatty acids (SCFAs), such as acetate, propionate, and butyrate. <i>Roseburia</i> can also produce butyrate from acetate promoting balance in energy homeostasis. The genus is widely recognized to influence colonic motility, support immunity, and suppress inflammation. Low levels are associated with several disease (including irritable bowel syndrome, obesity, Type 2 diabetes, nervous system conditions and allergies).



Therapeutic Options for Abnormally Low Commensal Bacterial Findings

- » Use a broad-spectrum, diverse probiotic formula, 50–450 billion
 CFUs/day depending on findings.
 May contain: Lactobacillus acidophilus,
 Bifidobacterium bifidum, Bifidobacterium longum, Lactobacillus rhamnosus,
 Bifidobacterium breve, Lactobacillus casei, Streptococcus thermophilus.
- » Increase dietary intake of vegetables and fibers (psyllium, oat bran)
- » Remove dietary sugar and refined carbohydrates
- » Prebiotic supplementation (resistant starch, xylooligosaccharide, inulin, beta-glucan, arabinogalactan)
- » Fermented foods, if tolerated
- » Reduce inflammation and address other imbalances on the GI-MAP

Therapeutic Options for Abnormally High Commensal Bacterial Findings

- » Consider any additional findings on GI-MAP and treat accordingly
- » Re-establish commensal bacteria using 5R protocol (see Table 2)
- » Remove dietary sugar and refined carbohydrates
- » In certain situations, overgrowth of commensal bacteria may be treated judiciously with antimicrobial herbs when all other findings are normal.

FIRMICUTES AND BACTEROIDETES PHYLA

Gram-negative *Bacteroidetes* and gram-positive *Firmicutes* are bacterial phyla that dominate the entire human digestive tract, including the mouth, nose, throat, and colon. An abnormal result in one or both of these phylum suggest imbalanced normal microbes in the GI tract. Further, high *Firmicutes* and low *Bacteroidetes* (*resulting in a high F/B ratio*) suggest microbial imbalance which may be related to increased caloric extraction from food, fat deposition and lipogenesis, impaired insulin sensitivity, and increased inflammation.

High Firmicutes/Bacteroidetes Ratio

- · Causes:
 - » Poor diet
 - » Dysbiosis
 - » Maldigestion or hypochlorhydria

- » Balance commensal bacteria using the 5R Protocol (see Table 2)
- » When Firmicutes phyla is high, consider using Bifidobacteria probiotics and Saccharomyces boulardii primarily.
- » Lactobacillus spp. and Bacillus spp. (found in probiotics) can elevate Firmicutes
- » Optimize the diet; a lower fat diet may help to normalize the F/B ratio
- » Address all other imbalances on the GI-MAP



OPPORTUNISTIC BACTERIA

Many bacteria measured on the GI-MAP are considered opportunistic pathogens, as they only cause disease and illness in some individuals, particularly the immune-compromised. Many individuals come into contact with opportunistic bacteria and experience no symptoms. Most sources consider these microbes to be normal in the stool. However, they can cause gastroenteritis and inflammation at high levels in vulnerable patients. Symptoms may include diarrhea, loose stools, abdominal pain, or even constipation. Overgrowth and excessive colonization by opportunistic bacteria may occur when the commensal bacteria are impaired by poor diet, antibiotic use, parasitic infection, or a weakened immune system. When intestinal permeability is present (see zonulin), these microbes could escape the lumen of the gut and infect extraintestinal sites.



OPPORTUNISTIC/OVERGROWTH MICROBES

Table 6.

DYSBIOTIC & OVERGROWTH BACTERIA		
Bacillus spp.	Common group of gram-positive bacteria in the <i>Firmicutes</i> phylum. Some strains are used as probiotics. High levels may result from reduced digestive function, SIBO, or constipation.	
Enterococcus faecalis Enterococcus faecium	Gram-positive species in the <i>Firmicutes</i> phylum. High levels may result from reduced stomach acid, PPI use, compromised digestive function, SIBO or constipation. High natural resistance to some antibiotics, which may result in overgrowth.	
Morganella spp.	Gram-negative group in the <i>Proteobacteria</i> phylum. May produce histamine. High levels may indicate increased intestinal inflammatory activity. High levels may cause diarrhea, and may also be associated with SIBO.	
Pseudomonas spp. Pseudomonas aeroginosa	Gram-negative bacteria in the <i>Proteobacteria</i> phylum. High levels may indicate increased intestinal inflammatory activity and may cause abdominal cramping and loose stools. Some strains of <i>P. aeroginosa</i> may produce toxins that can damage cells.	
Staphylococcus spp. Staphylococcus aureus	Gram-positive bacteria in the <i>Firmicutes</i> phylum. High levels may result from reduced digestive capacity, and intestinal inflammatory activity. Some strains may produce toxins and contribute to loose stools or diarrhea.	
Streptococcus spp.	Gram-positive bacteria in the <i>Firmicutes</i> phylum. <i>Streptococcus</i> spp. colonize skin and mucous membranes throughout the body; High levels in the intestine may result from low stomach acid, PPI use, reduced digestive capacity, SIBO or constipation; Elevated levels may also be indicative of intestinal inflammatory activity, and may cause loose stools.	

COMMENSAL & OVERGROWTH MICROBES		
Desulfovibrio spp.	A genus of Gram-negative sulfate reducing bacteria. The bacteria produce hydrogen sulfide (H2S), a metabolite which can influence cell signaling and reduce oxidative stress at low concentrations and pose toxicity at higher concentrations.	
Methanobacteriaceae (family)	Family of bacteria-like microbes that produce methane. Facilitates carbohydrate fermentation and short-chain fatty acid production by beneficial bacteria. High levels linked to chronic constipation, as well as some types of SIBO and IBS. Low levels may indicate reduced production of short-chain fatty acids and may be associated with inflammation.	

INFLAMMATORY & AUTOIMMUNE-RELATED BACTERIA		
Citrobacter spp. Citrobacter freundii.	Gram-negative bacteria in the <i>Proteobacteria</i> phylum. High levels may indicate increased intestinal inflammatory activity.	
Klebsiella spp. Klebsiella pneumoniae	Gram-negative bacteria in the <i>Proteobacteria</i> phylum. Common residents of the oral cavity and respiratory tract. May cause diarrhea, gas, abdominal pain, and bloating; Common after long-term antibiotic use; May release histamine in the gut; High levels may indicate increased intestinal inflammatory activity.	
Mycobacterium avium subsp. paratuberculosis	Bacterial species in the <i>Actinobacteria</i> phylum. Higher levels have been associated with Crohn's disease and rheumatoid arthritis.	
Proteus spp. Proteus mirabilis	Gram-negative bacteria in the <i>Proteobacteria</i> phylum. High levels may indicate increased intestinal inflammatory activity; May contribute to loose stools or diarrhea; Pets or wild animals can be a source.	

Table 6a.

COMMENSAL INFLAMMATORY & AUTOIMMUNE-RELATED BACTERIA		
Enterobacter spp.	Gram-negative genus in the <i>Proteobacteria</i> phylum. Closely related to <i>E. coli (in the same taxonomic family)</i> . High levels may indicate increased intestinal inflammatory activity. Low levels may indicate reduced mucosal health.	
Escherichia spp.	Gram-negative genus in the <i>Proteobacteria</i> phylum. Normal gut flora. <i>Escherichia coli</i> (<i>E. coli</i>) is the primary species in this genus. Most <i>E. coli</i> are nonpathogenic (pathogenic <i>E. coli strains are measured separately in "Pathogens" section of the GI-MAP).</i> High levels may be indicative of increased intestinal inflammatory activity. Low levels may indicate reduced mucosal health and decreased protection against pathogenic <i>E. coli</i> .	
Fusobacterium spp.	Genus of gram-negative bacteria in the <i>Fusobacteria</i> phylum. Commonly found in the oral cavity, and may also be found in the intestine. Associated with inflammatory processes, as well as autoimmune conditions such as systemic sclerosis.	
Prevotella spp.	Gram-negative species in the <i>Bacteroidetes</i> phylum. Associated with rheumatoid arthritis. High levels may result from reduced digestive capacity, or a high-starch diet.	

- Therapeutic Options and Considerations for Abnormally High Levels of Opportunistic Bacteria
 - Consider high-dose probiotics (300+ billion CFU/d)
 - Consider broad-spectrum antimicrobial herbs including: berberine, caprylic acid, garlic oil, oil of oregano, uva ursi, or olive leaf extract
 - » Optimize diet (low sugar, low refined carbs, high plant-based foods and fiber)
 - » See SIgA level to determine mucosal immunity and if patient is protected from overgrowth symptoms
 - » Use the 5R Protocol (see Table 2)
 - » Identify and remove potential sources of contamination or re-infection
 - » Address all other imbalances on the GI-MAP

- » If using antibiotics, see the Physician's Desk Reference for appropriate antibiotics for the specific microorganisms that are overgrown
- » If using antibiotics, consider rifaxamin, which remains in the GI tract and is also used to treat small intestinal bacterial overgrowth (SIBO)



Bacillus cereus



OPPORTUNISTIC BACTERIA

- Opportunistic Bacteria as a Trigger for Autoimmunity
 - » Certain opportunistic bacteria may initiate autoimmune thyroiditis or inflammatory arthritis such as rheumatoid arthritis and ankylosing spondylitis. These bacteria may trigger

or sustain the autoimmune process. Gastrointestinal symptoms are less common when these bacteria are elevated. When intestinal permeability is present (see zonulin), these microbes could escape the lumen of the gut and infect extraintestinal sites.

Table 7. Opportunistic Bacteria and Viruses Associated with Autoimmunity.

OPPORTUNISTIC BACTERIA AND AUTOIMMUNE ASSOCIATION		
Citrobacter spp.	Rheumatoid arthritis	
Fusobacterium spp.	Systemic sclerosis or inflammatory bowel disease	
Klebsiella spp.	Crohn's disease, ulcerative colitis, ankylosing spondylitis, and other spondyloarthropathies (which include ankylosing spondylitis, arthritis associated with Crohn's or ulcerative colitis, psoriatic arthritis, and reactive arthritis)	
M. avium subsp. paratuberculosis	Rheumatoid arthritis, Crohn's disease, Type I diabetes, possibly psoriasis	
Prevotella copri	Rheumatoid arthritis	
Proteus spp.	Rheumatoid arthritis	
Proteus mirabilis	Rheumatoid arthritis and spondyloarthropathies (listed above)	

VIRUSES AND AUTOIMMUNE ASSOCIATION		
Cytomegalovirus (CMV)	Systemic lupus erythematosus, systemic sclerosis, type 1 diabetes, rheumatoid arthritis	
Epstein Barr Virus (EBV)	Rheumatoid arthritis, lupus, Sjogren's, multiple sclerosis, autoimmune thyroid disorders	




FUNGI/YEAST

Fungal organisms are commonly found in the human digestive tract, but fungal overgrowth can cause illness in susceptible individuals. Fungal growth may be localized in the body. For instance, *Candida* spp. may be high in the large intestine but normal in the small intestine, and vice versa. In a patient with suspected fungal overgrowth, additional tests may be necessary to understand the complete picture of fungal overgrowth. Urinary D-arabinitol or antibodies to *Candida* are sometimes used.

- Common Causes of Yeast Overgrowth Include:
 - » Antibiotic use
 - » High intake of sugar, starches, and dietary fungi (beer, bread, nuts, cheese, corn)
 - » Hypochlorhydria
 - » Impaired immune function
 - » Dysbiosis
- Fungi/Yeast Targeted on the GI-MAP
 - » Candida albicans and Candida spp.
 - Commensal fungi that
 can be pathogenic to
 immunocompromised patients.
 Causes vaginal yeast infections
 and can be fatal in systemic
 infections. May cause diarrhea.

- Has been suggested to cause a cluster of symptoms including GI complaints, fatigue, and muscle or joint pain but evidence is weak.
- » Geotrichum spp.
 - May cause disease in immunosuppressed patients.
 Low levels may be a dietary artefact; certain strains are used to make soft cheeses.
- » Microsporidia spp.
 - The GI-MAP specifically detects
 Encephalitozoon intestinalis, which affects the GI. May cause diarrhea and wasting. Can disseminate to ocular, genitourinary, and respiratory tracts.





Candida albicans

- » Rhodotorula spp.
 - Common in soil, plants, bathrooms, and in beverages like milk, juice, and water. May be a commensal. Can cause disease in immunosuppressed patients.
- Common Symptoms of Fungal Dysbiosis
 - » GI symptoms: Gas, bloating, constipation, nausea, vomiting, and diarrhea.
 - » Other symptoms: Eczema, athlete's foot, vaginal yeast infections, thrush, and jock itch.

- Therapeutic Options and Considerations for Abnormally High Levels of Fungi/Yeast
 - » Reduce intake of sugars, starches, and fungi
 - » See SIgA levels and consider immune support
 - » Consider high-dose probiotics, Saccharomyces boulardii, and the 5R Protocol (see Table 2)
 - » Consider antifungal herbs such as caprylic acid, undecylenic acid, oregano oil, berberine, and/or garlic
 - » Consider pharmaceutical antifungals in severe cases. Nystatin is preferred because it stays in the GI tract.





VIRUSES

Cytomegalovirus

Epidemiology

- » Herpes virus that has infected 60% of the US population
- » One in three children have contracted CMV by five years old
- » Passed around in child daycare centers

Clinical Implications

- » Positive CMV on the GI-MAP indicates active infection of the GI, NOT past infection
- » Active infection may be asymptomatic or cause mild flu-like symptoms
- » CMV can also cause viral pneumonia, transaminitis, splenomegaly, colitis, fever, and encephalitis
- » Common in inflammatory bowel disease, immunocompromised patients
- » CMV colitis has a similar presentation to Clostridium difficile infection
- » CMV has been implicated in autoimmune diseases: lupus, systemic sclerosis, type 1 diabetes, and rheumatoid arthritis



Cytomegalovirus is a viral genus of the viral family known as Herpesviridae

Therapeutic Options and Considerations

- » No treatment is needed if asymptomatic
- » Prevent spreading CMV with regular handwashing
- » Antiviral herbs such as cat's claw, osha root, reishi mushrooms, vitamins A, C, and D, zinc, Echinacea
- » Address other imbalances on the GI-MAP and use 5R Protocol (see Table 2) to rebuild gut health and gut immunity



Epstein Barr Virus

- Epidemiology
 - One of the most common viruses worldwide; infects
 90–95% of the population
 - » Commonly contracted in childhood and causes mild symptoms
- Clinical Implications
 - » Positive finding on the GI-MAP indicates active EBV infection of the GI, not past infections
 - » Can cause infectious mononucleosis (mono)
 - » Symptoms include fatigue, fever, swollen lymph nodes, inflamed throat, enlarged spleen, and more
 - » May last two to four weeks in adolescents and adults
 - » May cause fatigue for weeks or months
 - » Associated with autoimmune conditions such as rheumatoid

- arthritis, lupus, Sjogren's, multiple sclerosis, autoimmune thyroid disorders
- » EBV may increase the risk of gastric cancer; especially if *H. pylori* present
- » May cause colitis
- » Found in 30-64% of IBD patients
- Therapeutic Options and Considerations
 - » Rest and hydration
 - » Antiviral herbs such as cat's claw, osha root; antiviral fungi such as reishi and/or Cordyceps mushrooms
 - » Vitamins A, C, and D, zinc, Echinacea
 - » Address other imbalances on the GI-MAP and use 5R Protocol (see Table 2) to rebuild gut health and gut immunity
 - » Follow-up blood testing may be indicated, including an EBV Early Antigen and EBV PCR test





PARASITES

A parasite is an organism that lives and feeds on a host organism at the expense of the host. The GI-MAP tests for pathogenic parasites and protozoa (some of which are non-pathogenic) most commonly occurring in the GI tract. Sources of exposure should be identified and eliminated to prevent reinfection.

PROTOZOA

Blastocystis hominis

- Epidemiology
 - » Fecal contamination of food or water is common
 - » Found worldwide
- Clinical Implications
 - » Symptoms include diarrhea, abdominal pain, nausea and vomiting, fever, fatigue, irritable bowel syndrome, infective arthritis
- Therapeutic options and considerations
 - » Difficult to eradicate

- » Consider nitazoxanide or tinidazole, oregano oil, and *S. boulardii*
- » Herbal treatments may not be as effective
- » Consider Artemisia, Coptis, or other broad-spectrum antiparasitic herbal formulas
- » Infection can be treated with metronidazole, iodoquinol or trimethoprim/sulfamethoxazole*
- » Consider probiotics and 5R Protocol (see Table 2)



Chilomastix mesnili

- Epidemiology
 - » Fecal contamination of food or water
- Clinical Implications
 - » Considered non-pathogenic and may not cause symptoms
 - » May indicate dysbiosis or suppressed immunity
- Therapeutic Options and Considerations
 - » Look for and address sources of fecal-oral contamination
 - » Consider probiotics and 5R Protocol (see Table 2)
 - » Address other imbalances on the GI-MAP

Cyclospora spp. (Cyclospora cayetanensis)

- Epidemiology
 - » Fecal contamination of food and water
 - » Associated with water- and food-borne outbreaks
 - » Common cause of traveller's diarrhea
 - » May be found on imported fresh produce from tropical regions
- Clinical Implications
 - Symptoms include prolonged watery diarrhea, abdominal cramping, loss of appetite, weight loss, nausea, and vomiting
 - » May cause alternating diarrhea and constipation



Cyclospora

- » Can cause bloating, flatulence, and burping
- » Flu-like symptoms such as fatigue, headaches, and low fever may be present in some individuals
- » Infection is usually self-limiting, with symptoms usually lasting about seven days, but can last weeks or months in immunosuppressed patients

Therapeutic Options and Considerations

- » In cases lasting more than seven days, treatment with an antibiotic combination of trimethoprim and sulfamethoxazole may be necessary*
- » Consider probiotics, broad-spectrum anti-parasitic herbal formula, and the 5R Protocol (see Table 2)
- » Look for and address sources of reinfection

Dientamoeba fragilis

- Epidemiology
 - » Not well understood; probably fecal contamination of food or water





Clinical Implications

- » May be asymptomatic
- » May cause diarrhea, abdominal pain, nausea, fever, fatigue, weight loss, appetite loss, and/or fatigue

Therapeutic Options and Considerations

- » "Moderate" amounts of DNA, that are not above the laboratory reference range, may cause symptoms and warrant treatment
- » Infection can be treated with iodoquinol, paromomycin, or metronidazole*
- » Consider probiotics, broad-spectrum anti-parasitic herbal formula, and the 5R Protocol (see Table 2)
- » Look for and address sources of reinfection
- » Address other imbalances on the GI-MAP

Endolimax nana

- Epidemiology
 - » Fecal contamination of food or water

• Clinical Implications

- » Considered non-pathogenic; individuals may be asymptomatic
- » May be indicative of dysbiosis, conservative treatment may be indicated if clinical presentation is consistent with enteroparasitosis

Therapeutic Options and Considerations

» Consider probiotics and the 5R Protocol (see Table 2)

- » Look for and address sources of fecal contamination
- » Address other imbalances on the GI-MAP

Entameoba coli

- Epidemiology
 - » Fecal contamination of food or water
 - » Found in the large intestine, considered to be non-pathogenic

Clinical Implications

- » May be indicative of dysbiosis, conservative treatment may be indicated if clinical presentation is consistent with enteroparasitosis
- Therapeutic Options and Considerations
 - » Consider probiotics and the 5R Protocol (see Table 2)
 - » Look for and address sources of fecal contamination
 - » Address other imbalances on the GI-MAP

Pentatrichomonas hominis

- Epidemiology
 - » Fecal contamination of food or water
- Clinical Implications
 - » Considered harmless, a non-pathogen
 - » Infected individuals are usually asymptomatic
 - » May contribute to dysbiosis
 - » Also colonizes dogs, cats, and other animals



Therapeutic Options and Considerations

- » May be asymptomatic
- » In women with vaginosis, consider treatment to reduce chances of vaginal contamination or reinfection (find treatments for Trichomonas vaginalis elsewhere)
- » If treatment is needed, consider a broad-spectrum antiparasitic herbal formula
- » Consider probiotics and the 5R Protocol (see Table 2)
- » Look for and address sources of fecal contamination
- » Address other imbalances on the GI-MAP

WORMS

Ancylostoma duodenale and Necatur americanus (Hookworms)

Epidemiology

- Infection occurs via skin contact
 with soil contaminated with
 larvae or ingestion of larvae
- » Infected cats and dogs are a source of exposure
- » Prevalent in southern Europe, Northern Africa, India, Asia, Caribbean islands, South America, and small areas of the United States
- » Associated with poor sanitation, inadequate housing construction, and lack of access to medications

Clinical Implications

- » Early symptoms are itching and rash where the larvae penetrated the skin
- » Symptoms of heavy infestations include: abdominal pain, diarrhea, fatigue, weight loss, iron deficiency anemia (IDA), coughing, and loss of appetite
- » Infected individuals may also be asymptomatic

Therapeutic Options and Considerations

- » Heavy infections can be treated with albendazole or mebendazole*
- » Individuals presenting with IDA may need iron supplementation
- » Consider anti-parasitic herbal treatments, gut immunity support, and the 5R Protocol (see Table 2)
- » Look for and remove sources of reinfection

Ascaris lumbricoides (Roundworm)

Epidemiology

- » Fecal contamination of food or water
- » Common in international travellers and recent immigrants from Latin America and Asia

· Clinical Implications

- » Early symptoms include fever, coughing, wheezing, and dyspnea
- » Late symptoms include abdominal pain, nausea, vomiting, frequent throat clearing, dry cough, "tingling throat," appendicitis, pancreatitis, and obstruction



- » Can cause reactive arthritis
- Therapeutic Options and Considerations
 - » Infections may be treated with albendazole, mebendazole, or ivermectin*
 - » Consider anti-parasitic herbal treatments, gut immunity support, and the 5R Protocol (see Table 2)
 - » Look for and remove sources of reinfection

Trichuris trichiura (Whipworm)

- Epidemiology
 - » Fecal contamination of produce or person-to-person contact
 - Prevalent in Asia, Africa,
 South America, and rural
 southeastern United States
- Clinical Implications
 - » Most individuals are asymptomatic, however diarrhea with mucus and blood may occur in some infected individuals
- Therapeutic Options and Considerations
 - » Infections may be treated with albendazole, mebendazole, or ivermectin*
 - » Individuals presenting with IDA may need iron supplementation



Microscopic cross section of a whipworm (Thichuris trichiura)

- » Consider anti-parasitic herbal treatments, gut immunity support, and the 5R Protocol (see Table 2)
- » Look for and remove sources of reinfection

Taenia spp. (Tapeworm)

- Epidemiology
 - » Fecal contamination of undercooked pork (*T. solium*) or beef (*T. saginata*)
 - » T. solium is found worldwide, but prevalent in communities who raise and eat pigs
 - » T. saginata is prevalent in Africa, parts of Eastern Europe, the Philippines, and Latin America where people raise cattle and eat raw beef

Clinical Implications

- » May be asymptomatic or present with mild symptoms
- » Symptoms include abdominal pain, nausea, weakness, increased appetite, loss of appetite, headache, constipation, dizziness, diarrhea, pruritus ani, hyperexcitability, and anemia

Therapeutic Options and Considerations

- » Infections may be treated with albendazole or praziquantel*
- » Consider anti-parasitic herbal treatments, gut immunity support, and the 5R Protocol (see Table 2)
- » Look for and remove sources of reinfection



INTESTINAL HEALTH MARKERS

DIGESTION

Pancreatic Elastase 1

Elastase 1 is a digestive enzyme secreted exclusively by the pancreas, giving a direct indication of pancreatic function. Elastase 1 is unaffected by pancreatic enzyme replacement therapy.

- Causes of Low Elastase 1:
 - » Suppressed pancreatic function
 - » Gallstones
 - » Hypochlorhydria, especially if *H. pylori* present
 - » Cystic fibrosis
 - » Low levels may be found in vegetarians/vegans
- Common Approaches for Addressing Low Pancreatic Digestive Enzyme Levels:
 - » Digestive support with betaine HCL

- » Chew thoroughly and relax at meal time
- » Pepsin
- » Plant or pancreatic enzyme supplements
- » Digestive herbs
- » Bile salts
- » Taurine
- » Consider underlying causes

Table 8. Staging of pancreatic insufficiency based on fecal elastase-1.

Fecal Elastase-1 Result	Clinical Significance
< 200 ug/g	Pancreatic insufficiency
200-500 ug/g	Decreased pancreatic output
> 500 ug/g	Normal pancreatic output





Steatocrit

Fecal fats are normally emulsified by bile salts and absorbed in the small intestines. High levels of fat in the stool may be an indication of maldigestion, malabsorption, or steatorrhea.

Causes of Elevated Steatocrit:

- » Hypochlorhydria
- » Maldigestion
- » Malabsorption
- » Pancreatic insufficiency (see elastase-1)
- » Bile salt insufficiency
- » Improper mastication
- » Celiac disease

Therapeutic Approaches and Considerations for High Fecal Fats:

- » Support digestion with betaine HCL
- » Pepsin
- » Digestive herbs or "bitters"
- » Bile salts
- » Taurine
- » Consider underlying causes of malabsorption, such as celiac disease, dysbiosis, or food sensitivities

ADDITIONAL GI MARKERS

Beta-Glucuronidase

High levels of fecal beta-glucuronidase can indicate unfavorable metabolic changes in the colon. Beta-glucuronidase may indicate dysbiosis and interference with Phase II detoxification involving glucuronidation.

• Major Producers of β-glucuronidase are:

- » Bacteroides fragilis
- » Bacteroides vulgatus
- » Bacteroides uniformis
- » Clostridium paraputrificum
- » Clostridium clostridioforme
- » Clostridium perfringens
- » Escherichia coli
- » Eubacterium
- » Peptostreptococcus
- » Ruminococcus
- » Staphylococcus

Clinical Indications of High β-glucuronidase:

- » Dysbiosis in the colon or small intestinal bacterial overgrowth (SIBO)
- » Extremely elevated cases associated with colon cancer risk
- » Problems with detoxification, especially estrogen (via glucuronidation pathway)
- » Overexposure to toxins or drugs

Therapeutic Approaches and Considerations for Elevated β-glucuronidase:

- » Address dysbiosis, if present
- » Promote bacterial diversity with probiotics, fiber, prebiotics, and fermented foods
- Consider liver support such
 as milk thistle and calcium
 D-glucarate, especially if patient
 is taking hormone replacement
 or has increased cancer risk
- » If there are no signs of dysbiosis on the GI-MAP, consider a SIBO breath test



Occult Blood Fecal Immunochemical Testing (FIT)

FIT is quantitative and directly measures the concentration of hemoglobin present in stool, rather than just the qualitative presence of hemoglobin. This test uses antibodies specific for human hemoglobin and therefore does not require dietary restrictions or multiple samples, significantly reducing the appearance of false positives. This method has better detection of lower hemoglobin concentrations than qualitative tests, eliminating potential false negatives as well. Literature suggests a result of 10 ug/g may be indicative of potentially more serious conditions such as polyps or colorectal cancer. A variety of ailments can cause lower counts of blood in stool, such as hemorrhoids, anal fissures, pathogenic infection such as giardia, liver disease, and upper GI infections.

Possible Causes of Positive Occult Blood:

- » Bleeding ulcer
- » Inflammatory bowel disease
- » Cancer
- » Intestinal polyps
- » Upper GI bleeds that cause iron deficiency anemia
- Common Approaches for Addressing Fecal Occult Blood
 - » Identify source
 - » Follow-up testing recommended

IMMUNE RESPONSE

Secretory IgA (SIgA)

Immunoglobulin A is the primary immunoglobulin in the intestinal mucosa. It represents a "first line of defense" in response to antigens and pathogens in the GI and respiratory tracts. In addition to protecting against pathogens, SIgA plays a major role in helping to maintain balance in the microbiome and protecting against exposure to food-derived antigens.

Low Fecal SigA – The gut immune system is suppressed. Investigate underlying causes, such as chronic dysbiosis, antigen exposure, chronic stress, immunocompromised patient, or even protein malnutrition.

- Therapeutic Approaches for Low SIgA Levels:
 - » Address any chronic GI infections, if appropriate
 - » Address microbiome imbalances
 - » Address chronic stress and adrenal health, if needed
 - » Colostrum or immunoglobulins
 - » Supplement with S. Boulardii
 - » GI mucosal support with glutamine
 - » Lactobacillus and Bifidobacteria probiotics
 - » General immune support
 - » Essential fatty acids
 - » Zinc
 - » Address other imbalances on the GI-MAP





High Fecal SIgA – Elevated immune response to antigens in the GI tract. Investigate underlying causes, such as chronic dysbiosis, acute infections, acute stress, or food sensitivities.

- Therapeutic Approaches for High SIgA Levels:
 - » Address GI infections
 - » Address any food allergies and sensitivities
 - » General immune support

Anti-gliadin SIgA – Gliadin is a component of gluten, the protein found in wheat and other field grass grains such as barley, malt, and rye. The presence of fecal anti-gliadin antibodies can indicate an immune response (in the gut) to gluten in the diet. Fecal anti-gliadin antibodies do not necessarily correlate with blood levels.

High Anti-gliadin SIgA – Elevated immune response to gliadin in the lumen of the gut.

- Treatment:
 - » Consider gluten elimination for a trial period
 - » If patients have been gluten-free, consider hidden sources of gluten and gliadin cross-reactive food such as dairy, corn, oats, millet, rice and yeast
 - » Consider intestinal barrier support, including supplements such as L-glutamine, zinc carnosine, and colostrum

Eosinophil Activation Protein (EDN/EPX)

EDN/EPX is a protein released by activated eosinophils which has strong cytotoxic characteristics. The protein plays a significant role in a variety of inflammatory and mast-cell mediated pathologies in addition to fighting pathogens, particularly viral infections.

High Eosinophil Activation Protein -

The accumulation of eosinophil activation protein in the intestine is associated with inflammation and tissue damage, and the level of EPA in the stool can serve as an objective measure for chronic inflammation in the GI tract. In the case of inflammatory bowel disease, the marker can be used to evaluate disease activity and predict relapse. The EPA marker can also be used to determine the effectiveness of a food elimination diet to control symptoms or disease progression.

- Possible Causes of Elevated Eosinophil Activation Protein:
 - » Respiratory allergies
 - » Asthma
 - » Food allergies and sensitivities
 - » IBD
 - » IBS
 - » Eosinophilic esophagitis (EE)
 - » Functional dyspepsia
 - » Acid reflux
 - » Intestinal barrier damage/dysfunction
 - » Anxiety (IBS-related anxiety)
 - » Intestinal parasites



INFLAMMATION

Calprotectin

Fecal calprotectin is the most studied marker of gastrointestinal inflammation. High calprotectin indicates neutrophil infiltration to the gut mucosa. Calprotectin is the gold standard marker for the diagnosis and monitoring of inflammatory bowel disease. It is used to differentiate IBD from irritable bowel syndrome.

- Possible Causes of Elevated Calprotectin
 - » Intestinal infections and proinflammatory dysbiosis
 - » Food allergens, toxins and certain drugs (e.g., non-steroidal antiinflammatory drugs [NSAIDs])
 - » Inflammatory bowel disease
 - » Polyps
 - » Diverticulitis
 - » Colorectal cancer
- Therapeutic Approaches and Considerations
 - » Address possible causes of elevated calprotectin
 - » Persistently elevated calprotectin may indicate chronic inflammatory disease; further evaluation by a qualified medical professional is advised
 - » Consider anti-inflammatory support (e.g., anti-inflammatory diet, curcumin, omega-3 fatty acids, aloe, and resveratrol)

ADD-ON TESTS

The following optional add-on tests can be ordered with the GI-MAP.

Zonulin

Zonulin is a protein that opens intercellular tight junctions in the gut lining (the connections between epithelial cells that make up the gastrointestinal lining). Zonulin increases intestinal permeability in the jejunum and ileum and is considered a biomarker for barrier permeability.

- Therapeutic Approach for Elevated Intestinal Permeability:
 - » Address dysbiosis (pathogens and opportunistic microbe overgrowth, lack of beneficial microbes)
 - » Eliminate gluten, address potential food sensitivities
 - » Promote a healthy intestinal barrier with L-glutamine, butyrate, essential fatty acids, aloe vera, probiotics, zinc carnosine, slippery elm, marshmallow, deglycyrrhizinated licorice

Zonulin is available as an optional add-on to the GI-MAP, or as a stand-alone test.



Gluten Peptide

The Gluten Peptide test detects the presence of gluten in stool. Gluten is broken into gliadin and glutenin subunits during the digestive process. A specific section of the gliadin protein, referred to as the 33-mer peptide,* is resistant to digestion in the human gastrointestinal (GI) tract but can be accurately detected in stool several days after the consumption of gluten.

The Gluten Peptide test can help practitioners monitor a patient's gluten-free (GF) diet compliance. Any detected level of the fecal gluten protein indicates dietary exposure to gluten within the previous 2–4 days.[†]

Clinical Implications:

- » For patients following a GF diet, detected levels of the gluten peptide can indicate accidental exposure from hidden sources or food cross-contamination
- For patients regularly consuming gluten, detected levels are expected findings

The Gluten Peptide marker is available as an optional add-on to the GI-MAP, or as a standalone test.

Universal Antibiotic Resistance (AR) Genes Panel

The Universal Antibiotic Resistance Genes Panel detects the presence of 55 genetic elements associated with resistance to 10 different classes of antibiotics. Useful for patients who have been hospitalized, treated with antibiotics, or who have stubborn, chronic infections, the Universal Antibotic Resistance Genes panel can inform practitioners if antibiotic resistance is a potential stumbling block for patients.

The Universal Antibiotic Resistance (AR) Gene Panel is available as an optional add-on only to the GI-MAP, or Pathogens Panel. See page 42 for more information regarding AR Genes testing on the GI-MAP and antibiotic stewardship.

StoolOMX™

StoolOMX measures gut microbial metabolites that evaluate 25 bile acids and 9 short chain fatty acids, providing actionable insights into conditions like IBS, IBD, and bile acid diarrhea.

StoolOMX is available as an optional addon to the GI-MAP. Please see StoolOMX on page 44 for a detailed interpretive overview.



^{*} Because of structural homology (the structural similarity) of gluten proteins, detected 33-mer gliadin protein may indicate the consumption of wheat, rye, barley, corn, or oat grain.

[†] Speed of digestion (motility and transit time) can influence the transience of gluten peptide in stool.

| ANTIBIOTIC | RESISTANCE | GENES

In an effort to promote antibiotic stewardship, the GI-MAP **includes** results of *H. pylori Antibiotic Resistance Genes* in the microbiome as part of the test. Additionally, the *Universal Antibiotic (AR) Resistance Genes* panel can be ordered as an **add-on** to the GI-MAP or to the GI-Pathogens panel.

Antibiotic Stewardship

Antibiotic stewardship is the effort to measure and improve how antibiotics are prescribed by clinicians and used by patients. If an antibiotic resistance gene is present, then that class of antibiotics is designated POSITIVE for antibiotic resistance. A positive result for the presence of resistance genes for a given antibiotic indicates that the antibiotic is not an ideal choice for an antibiotic protocol.

H. pylori Antibiotic Resistance Genes

The *H. pylori* Antibiotic Resistance Genes are specific to the *H. pylori* genome detected on an individuals' GI-MAP. These results will only be reported if the *H. pylori* result is 5.0e2 or greater.

Universal Antibiotic Resistance Genes

The Universal Antibiotic resistance genes apply to *all* of the microorganisms within the fecal sample. Since microbes can rapidly share DNA under stress, the presence of antibiotic resistance in any organism is reason enough to avoid that drug class.





	Result	Reference		Result	Reference
Amoxicillin	Positive	Negative	Clarithromycin	Positive	Negative
Genes associated with amoxic	illin resistance		Genes associated with clarithron	nycin resistance	
PBP1A S414R	Present		A2142C	Absent	
PBP1A T556S	Absent		A2142G	Absent	
PBP1A N562Y	Absent		A2143G	Present	
	Result	Reference		Result	Reference
Fluoroquinolones	Negative	Negative	Tetracycline	Negative	Negative
Genes associated with fluoroqu	uinolone resistance		Genes associated with tetracyclii	ne resistance	
gyrA N87K	Absent		A926G	Absent	
gyrA D91N	Absent		AGA926-928TTC	Absent	
gyrA D91G	Absent				
gyrB S479N	Absent				
gyrB R484K	Absent				

Figure 4. Sample *H. pylori* Antibiotic Resistance Genes Section. *See additional details in the "Helicobacter pylori" section of this guide (Page 17, Figure 2).*

		UNIVER	SAL ANTIBIO	TIC RES	STANCE	GENES		
b-Lactams	Result Negative	Reference Negative	Macrolides	Result Negative	Reference Negative	Trimethoprim	Result Negative	Reference
D Luciumo	Hogunie	Hogalive	Muoronaco	Hogunie	Hogunie	Timethopimi	Hogunie	Negative
blaNDM-1	Absent		acrA	Absent		dfrA1	Absent	
CTX-M 1	Absent		acrB	Absent		dfrA12	Absent	
CTX-M 2	Absent		emrE	Absent		dfrA14	Absent	
CTX-M 8/25	Absent		ermA	Absent		dfrA15	Absent	
CTX-M 9	Absent		ermB	Absent		dfrA17	Absent	
GES	Absent		ermC	Absent		dfrA5	Absent	
OXA-1	Absent		macA	Absent		dfrA7	Absent	
PER-1	Absent		macB	Absent		dfrB1	Absent	
PER-2	Absent		mefA	Absent		dfrB2	Absent	
SHV	Absent		mphA	Absent		dfrB3	Absent	
TEM	Absent		msrA	Absent				
VEB	Absent		tolc	Absent				
Fluoroquinolones	Positive	Negative	Ciprofloxacin	Negative	Negative	Sulfonamides	Positive	Negativ
gnrA	Present		emea	Absent		sul1	Present	
gnrB	Absent		pmra	Absent		sul2	Present	
qnrS1	Absent		·			sul3	Present	
Vancomycin	Negative	Negative	Nitroimidazoles	Negative	Negative	Methacillin	Negative	Negativ
vanA	Absent		nimA	Absent		mecA	Absent	
vanA2	Absent		nimB	Absent		Chloramphenicol	Negative	Negativ
vanB	Absent		nimC	Absent		_		J
vanC1	Absent		nimD	Absent		cata13	Absent	
vanC2-1	Absent		nimE	Absent				
vanC2-2	Absent							

Figure 5. Sample Universal Antibiotic Resistance Genes Panel. *Available as an optional add-on only to the GI-MAP, or Pathogens Panel.*



$STOOLOMX^{TM}$

StoolOMX measures gut microbial metabolites—a comprehensive set of bile acids (BA) and short chain fatty acids (SCFA) via LC-MS/MS. This data offers an in-depth understanding of the stool metabolome and its clinical applications. The StoolOMX profile can be ordered as an **add-on** to the GI-MAP.

BILE ACIDS

Primary bile acids are natural products produced in the liver from cholesterol synthesis. These primary bile acids are conjugated with either taurine or glycine to increase solubility and are stored in the gallbladder as bile.

Primary, conjugated bile acids are the main component in bile. At mealtime, they aid in the emulsification and absorption of dietary fats in the small intestine. After contributing to fat absorption, the majority (~95%) of primary bile acids are reabsorbed in the distal ileum and returned to the liver via the portal vein—a process called enterohepatic circulation. A small portion (~5%) of primary bile acids will reach the colon, where they are metabolized by gut bacteria and deconjugated (glycine and taurine removed) to produce unconjugated secondary bile acids. Deconjugation is generally favorable.

Commensal gut bacteria deconjugate bile acids to facilitate their reabsorption and recycling while also preparing them for further bacterial modification.

Under normal physiological conditions, most bile acids are reabsorbed, and very few bile acids are excreted in stool. Under abnormal physiological conditions, excess bile acids enter the colon and can promote inflammation and diarrheal symptoms while exerting negative effects on commensal bacteria.

Measuring concentrations and ratios of bile acids in stool can offer root-cause insights into digestive symptoms, malabsorptive disorders, immune system regulation, and even metabolic impacts. Furthermore, altered patterns of bile acid metabolites have emerging disease state associations and can be used as part of diagnostic workup and treatment management in conditions such as bile acid diarrhea (BAD), inflammatory bowel disease (IBD), and irritable bowel syndrome (IBS) (-C, -D, -M).

StoolOMX measures 25 bile acid metabolites in total concentrations, percents, and ratios. This data offers an in-depth understanding of the stool metabolome and impactful clinical considerations.



BILE ACIDS - SUMMARY

Total Bile Acids Concentration

- An absolute concentration of total bile acids in stool measured in ng/g
 - » Reflects the amount of bile acids reaching the colon and not reabsorbed in the distal ileum after digestion
- Elevated Total Bile Acids:
 - » Excess/elevated bile acids in stool is an indication of bile acid malabsorption (BAM) and often results in symptoms such as diarrhea (bile acid diarrhea – BAD)
 - » Quantification of bile acids may be part of a comprehensive workup for irritable bowel syndrome with diarrhea (IBS-D), as it is currently estimated that up to 1/3 of IBS-D patients present with BAD
 - » Elevated stool bile acids may suggest altered bile acid production or reabsorption, which could be due to liver disease, gallbladder dysfunction, intestinal issues, or diet composition

Low Total Bile Acids:

» Low total bile acid concentrations in stool are generally favorable but may be indicative of slow transit or constipation

Total Primary Bile Acids Percent

- Percent concentration of primary bile acids in stool
 - » Primary bile acids should be present in very low concentrations in stool
 - » Primary bile acids in the colon can be damaging to the gut lining and can have antimicrobial action on the microbiome
 - » Higher levels of primary bile acids in a stool sample are indicative of excess production or BAM and are associated with an array of GI symptoms and disease states, including IBD
 - See individual primary bile acid concentrations on page 2 of the StoolOMX report for further analysis and associations

Total Secondary Bile Acids Percent

- Inverse of total primary bile acids in stool
 - » Listed as percent concentration of secondary bile acids in stool
 - » Secondary bile acids should be present in high concentrations in stool
 - » Secondary bile acids play diverse roles in regulating intestinal motility, gut barrier function, metabolism, and immune balance
 - » In general, a higher number of secondary bile acids relate to a healthier, more diverse microbiome and normal physiological GI function
 - » See individual secondary bile acid concentrations on page 2 of the StoolOMX report for further analysis and associations



Table 9. This section represents the most abundant bile acids.

PERCENTAGE BREAKDOWN OF MOST ABUNDANT BILE ACIDS					
SECONDARY BILE ACIDS	FUNCTION	AVERAGE PERCENTAGE [†]			
Deoxycholic Acid (DCA)	DCA is a major secondary bile acid that is formed from the bacterial metabolism of CA. Low levels observed in ulcerative colitis (UC). Elevated levels are associated with liver and colon cancer and metabolic imbalance. Implicated in modulating the immune response by inhibiting pro-inflammatory cytokine production. Can inhibit <i>Lactobacillus</i> spp., <i>Bifidobacterium</i> spp., and <i>Bacteroides fragilis</i> .	48%			
Lithocholic Acid (LCA)*	LCA is a major secondary bile acid that is formed from the bacterial metabolism of CDCA. Low levels observed in UC. Elevated levels are associated with liver and colon cancer and may contribute to blood sugar dysregulation. Modulates immune responses at normal levels.	27%			
Isolithocholic Acid (Iso-LCA)	Prevalent secondary bile acid in stool that has positive associations with longevity. Influences metabolic health, interacts with cellular receptors, and modulates immune responses at normal levels. Influences the differentiation and function of T cells, particularly T helper 17 (Th17) cells, which play a critical role in the inflammatory response.	8%			
Other	The remaining primary and secondary bile acid metabolites measured on StoolOMX.	17%			

^{*} On StoolOMX, the LCA value is a summation of LCA + Allo-LCA.; † Reference set at 50th percentile.

LCA/DCA Ratio

LCA and DCA are secondary bile acids formed from CDCA and CA in the colon. CDCA is converted to LCA, and CA is converted to DCA. This ratio can be useful in determining the risk of certain disease states and conditions.

LCA is thought to be more toxic than DCA due to its inhibitory effects on antioxidant pathways. An elevated ratio can be seen in patients with gallstones and after cholecystectomy. A higher ratio is also associated with an increased risk for colon cancers, while a lower ratio may indicate a reduced cancer risk.

Bile Acids: Therapeutic Applications

Bile acids and the microbiome have a bidirectional relationship and impact. It is important to first evaluate patient symptoms, such as stool frequency and Bristol stool type.

Always analyze StoolOMX results with GI-MAP findings and support accordingly. It is important to assess liver function, as bile acids are originally produced in the liver.

If there are significant imbalances in the bile acids summary, review the list of individual bile acids on page 2 of the StoolOMX report.





- A healthy microbiome is typically associated with:
 - » Total bile acid concentration in normal range
 - » High secondary bile acid percent
- An unhealthy microbiome is typically associated with:
 - » Elevated total bile acid concentration
 - » High primary bile acid percent, as the microbiome is responsible for converting primary bile acids to secondary bile acids
- Common GI-MAP patterns that may be associated with abnormal StoolOMX results:
 - » Presence of pathogens
 - » Insufficiency dysbiosis
 - » Increased Firmicutes:Bacteroidetes ratio
 - » Overgrowth of GI-MAP (report page 3) inflammatory opportunists
 - » Elevated Steatocrit
 - » Decreased Elastase-1
 - » Elevated Calprotectin
- Lifestyle Considerations
 - Evaluate macronutrient composition of diet with respect to calories from fat.
 Patients may need to adopt a lower fat diet if StoolOMX results are abnormal.
 - » Evaluate and ensure adequate dietary fibers
 - » Anti-inflammatory dietary practices
 - » Support liver and gallbladder
 - » Incorporate regular exercise

Medication/Supplement Considerations for BAM

- » Bile acid sequestrants
 - ex: cholestyramine, colesevelam
- » Tauroursodeoxycholic acid (TUDCA)
- » Digestive support as needed
 - Digestive enzymes, betaine HCL, digestive bitters, ox bile, etc.
- » Activated charcoal
- » Pectin fibers
- » Probiotics
- » Polyphenols such as resveratrol

Table 10. GI-MAP microbes involved in deconjugation of primary bile acids to secondary bile acids.

GI-MAP MICROBES
Firmicutes phyla
Bacteroidetes phyla
Escherichia spp.
Bacteroides spp.
Bifidobacterium spp.
Enterococcus spp.
Lactobacillus spp.
Methanobacteriaceae family






Figure 6. Bile Acids and Fatty Acids Overview.

Primary bile acids are synthesized from cholesterol in the liver and conjugated with either taurine or glycine. They are stored in the gallbladder and released during digestion to assist with the absorption of fat and fat-soluble vitamins. Normally, ~95% of primary bile acids are reabsorbed via the portal vein, while ~5% are metabolized by gut bacteria to produce secondary bile acids.

Saccharolytic short chain fatty acids (SCFAs) are primarily metabolites of dietary fiber fermentation in the gut while proteolytic branched chain fatty acids (BCFAs) are metabolites of protein fermentation. Acetate, propionate, and butyrate are three major SCFAs, which account for ~90% of the SCFAs produced by gut microbiota. SCFAs are known to have numerous health effects and can enhance fecal excretion of bile acids.





SHORT CHAIN FATTY ACIDS – SUMMARY

Short chain fatty acids (SCFAs) are small (0-6 carbon) saturated fatty acids that are produced from the microbiome (gut bacteria) by fermentation of food substrates. Fermentation of carbohydrate—called saccharolytic fermentation—produces straight chain fatty acids (SCFA). Fermentation of protein—called proteolytic fermentation—produces branched chain fatty acids (BCFA).

StoolOMX measures 9 short chain fatty acid metabolites in percentages, total concentrations, and ratios. The test also includes the ratio of straight chain fatty acids to branched chain fatty acids. Ratios can provide insight into dietary composition and digestive impacts such as saccharolytic and proteolytic (putrefactive) fermentation and digestive insufficiency. Assessing SCFA levels in stool offers a deeper look into digestive health and the overall balance of the gut microbiome.

Major SCFAs

Acetate, butyrate, propionate, and valerate are the major primary straight chain metabolites. As preferred fuel for intestinal cells, these metabolites provide a variety of beneficial effects for intestinal health.

SCFAs serve as an energy source for colon cells, strengthen the gut barrier and provide metabolic and immune system signaling. They have a profound anti-inflammatory effect by inducing and selectively expanding T-regulatory cells (T regs) in the large intestine, which, in turn, suppresses the proinflammatory action of Th17. As metabolites of bacterial fermentation, their levels can reflect gut microbiota composition.

Low levels of SCFAs are linked to various conditions, including irritable bowel syndrome, obesity, and inflammatory bowel disease. Low levels may also be related to slow motility.

Table 11.

PERCENT BREAKDOWN OF MAJOR SHORT CHAIN FATTY ACIDS REPRESENTED ON STOOLOMX					
MAJOR SCFA	FUNCTION	AVERAGE PERCENTAGE†			
Acetate	Most abundant SCFA. Involved in lipid synthesis and appetite regulation, maintains energy balance and metabolic homeostasis. Resists oxidation and mitochondrial stress.	52%			
Butyrate	Primary energy source for colonic cells. Supports intestinal barrier integrity, reduces intestinal inflammation, promotes motility, enhances fatty acid oxidation, inhibits tumor cell progression, and fosters a balanced microbiome.	20%			
Propionate	Supports intestinal barrier integrity, impacts energy balance, gluconeogenesis, and lipid metabolism. Involved in appetite regulation. Proposed as a biomarker for IBS.	25%			
Valerate	Stimulates intestinal epithelium growth, inhibits colon cancer cell production, modulates glucose and lipid metabolism. Antimicrobial effects against <i>C. difficile</i> .	3%			

[†] Reference set at 50th percentile.



SCFA/BCFA Ratio

Saccharolytic straight chain fatty acids (SCFAs) are primary metabolites of dietary fiber fermentation in the gut while proteolytic branched chain fatty acids (BCFAs) are metabolites of protein fermentation.

Ideally, saccharolytic SCFAs should make up ~95% of total SCFAs, while BCFAs should make up ~5% of total SCFAs.

The production of branched chain fatty acids leads to other fermentation products that can be harmful to the colon epithelium, such as ammonia, phenol, p-cresol, or biogenic amines. BCFAs are mainly produced during fermentation of branched chain amino acids (valine, leucine, and isoleucine) by the intestinal microbiota.

High levels of BCFAs compared to saccharolytic SCFAs can indicate weak digestion, poor transit, and inflammatory dysbiosis. There may also be associations between elevated BCFA concentrations in stool and obesity, IBD, hypercholesterolemia, and metabolic-associated fatty liver disease (MAFLD).

The SCFA/BCFA ratio may decline with age, primarily due to a significant decrease in SCFA levels in stool.

Short Chain Fatty Acids: Therapeutic Applications

- Causes for Low SCFA Levels:
 - » Diarrhea (rapid transit leading to decreased SCFA production)
 - » Constipation (slow transit leading to increased SCFA absorption)
 - » Inflammation (high calprotectin)
 - » Chronic antibiotic use
 - » Low complex carbohydrate and fiber intake
 - » Insufficiency dysbiosis
 - » Inflammatory bowel disease
- Cause for High SCFA Levels:
 - » Increased transit time (diarrhea)
- · Causes for High BCFA Levels:
 - » High protein intake and low fiber and/or polyphenol intake (Western diet)
 - » Low total SCFAs
 - » Weak digestion/hypochlorhydria
 - » Increased age
 - » Metabolic imbalance
- Therapeutic Options and Considerations to Increase Levels of SCFAs:
 - » Ensure diet is high and diverse in plant-based fibers, polyphenols, and fermented foods
 - » Avoid antibiotics
 - » Consider probiotics to support butyrate-producing organisms
 - » Consider support with butyrate salts
 - » Intestinal barrier support
 - » Support digestion if indicated
 - » Maintain a healthy weight





Table 12. The bacteria involved in the production of saccharolytic straight chain fatty acids and proteolytic branched chain fatty acids. Acetate, propionate, and butyrate make up \sim 90–95% of fecal SCFAs, with valerate and caproate comprising a smaller percentage. Isobutyrate, isovalerate, 2-methylbutyrate, and isocaproate make up \sim 5–10% of fecal SCFAs. Microbiota listed below each metabolite are involved in the production of that respective fatty acid.

	BACTERIA INVOLVED IN SHORT CHAIN FATTY ACID PRODUCTION					
	SACCHAROLYTIC STRAIGHT CHAIN FATTY ACIDS (SCFA)	PROTEOLYTIC BRANCHED CHAIN FATTY ACIDS (BCFA)				
	Acetate • Blautia hydrogenotrophica • Bifidobacterium spp. • Lactobacillus spp. • Clostridium spp. • Streptococcus spp.	Isobutyrate • Bacteroides spp. • Clostridium spp.				
Percent of Fecal SCFAs	Propionate • Bacteroidetes • Negativicutes • Megasphaera elsdenii • Lachnospiraceae • Coprococcus catus • Akermansia municiphilia	Isovalerate • Bacteroides spp. • Clostridium spp.	t of Fecal SCFAs			
~90-95 Percer	Butyrate • Firmicutes phyla • Faecalibacterium prausnitzii • Roseburia spp. • Eubacterium spp. • Clostridium coccoides	2-Methylbutyrate (from leucine)Bacteroides spp.Clostridium spp.	~5-10 Percent of			
	Valerate (potentially toxic) • Clostridia					
	Caproate (potentially toxic - from lactate) • Megasphaera elsdenii • Clostridium spp. BS-1	Isocaproate (associated with disease)				

- Therapeutic Options and Considerations to Decrease Levels of SCFAs:
 - » Low-FODMAP diet
 - » Maintain a healthy weight
 - » Support digestion if indicated
 - » Consider multi-strain probiotics to decrease acetate concentrations
- Therapeutic Options and Considerations to Decrease Levels of BCFAs:
 - » Increase fiber intake (particularly insoluble fibers) and/or decrease protein intake
 - » Consider underlying associations such as inflammatory bowel disease, chronic kidney disease, obesity, type 2 diabetes mellitus, MAFLD, or impaired digestion



INDIVIDUAL BILE ACID RESULTS

Research is more widely available on the clinical implications of certain primary and secondary bile acids compared to others.

Table 13.

PRIMARY BILE ACIDS			
MARKER	ABBR	GUIDE	
Cholic Acid	CA	Major primary bile acid. May be increased in IBS-D patients. CA levels in stool may serve as a biomarker for IBS and UC. During a UC flare, deconjugation of CA to DCA may be impaired. Can inhibit <i>Roseburia</i> spp. and <i>Lactobacillus</i> spp.	
Chenodeoxycholic Acid	CDCA	Major primary bile acid. May be increased in IBS-D patients and may contribute to visceral hypersensitivity. Elevated levels may be associated with metabolic imbalance.	
Taurochenodeoxycholic Acid	TCDCA	Minor primary bile acid.	
Taurocholic Acid	TCA	A whole grain diet may increase levels of TCA.	
Glycochenodeoxycholic Acid	GCDCA	May be increased in IBS-D patients.	
Glycocholic Acid	GCA	A whole grain diet may increase levels of GCA.	
Hyocholic Acid	НСА	Lower levels are associated with pre-diabetes. HCA levels may serve as an indication of metabolic health.	

Table 14.

SECONDARY BILE ACIDS			
MARKER	ABBR	GUIDE	
Lithocholic Acid*	LCA	LCA is a major secondary bile acid that is formed from the bacterial metabolism of CDCA. Low levels observed in UC. Elevated levels are associated with liver and colon cancer and may contribute to blood sugar dysregulation. Modulates immune responses at normal levels.	
Deoxycholic Acid	DCA	DCA is a major secondary bile acid that is formed from the bacterial metabolism of CA. Low levels observed in ulcerative colitis (UC). Elevated levels are associated with liver and colon cancer and metabolic imbalance. Implicated in modulating the immune response by inhibiting pro-inflammatory cytokine production. Can inhibit <i>Lactobacillus</i> spp., <i>Bifidobacterium</i> spp., <i>Bacteroides fragilis</i> , and <i>Clostridium difficile</i> .	
Isolithocholic Acid	Iso-LCA	Prevalent secondary bile acid in stool that has positive associations with longevity. Influences metabolic health, interacts with cellular receptors, and modulates immune responses at normal levels. Influences the differentiation and function of T cells, particularly T helper 17 (Th17) cells, which play a critical role in the inflammatory response.	

 Table 14. Continued.

SECONDARY BILE ACIDS				
MARKER	ABBR	GUIDE		
12-Ketolithocholic Acid	12-KLCA	Prevalent secondary bile acid in stool. May provide anti-inflammatory effects in IBD, modulate cholesterol metabolism, and improve glucose homeostasis.		
3-oxoDeoxycholic Acid	3-oxoDCA	Prevalent secondary bile acid in stool that has positive associations with longevity.		
Ursodeoxycholic Acid	UDCA	Enhances bile flow, hepatoprotective, and has immunomodulatory properties. Levels in stool may serve as a biomarker of IBS and UC. Elevated levels may be associated with metabolic imbalance.		
Glycolithocholic Acid	GLCA	May be negatively correlated with fecal calprotectin levels in UC. Influences metabolic signaling pathways, modulates gut microbiome, involved in lipid absorption, and aids in the detoxification of potentially harmful bile acids.		
Glycoursodeoxycholic Acid	GUDCA	Modulates the gut microbiome, cytoprotective, and influences lipid and glucose metabolism. May be elevated in IBS-D.		
Glycodeoxycholic Acid	GDCA	Modulates the gut microbiome. Involved in lipid absorption and metabolic regulation. May be reduced in UC.		
Taurolithocholic Acid	TLCA	Involved in lipid absorption and metabolic regulation. Modulates inflammation and may exert cholestatic effects when dysregulated. A whole grain diet may increase levels of TLCA.		
Tauroursodeoxycholic Acid	TUDCA	Often used synthetically as a supplement. Used clinically for cholestatic liver diseases. Cytoprotective, particularly in the liver. Involved in the regulation of bile acid metabolism and lipid metabolism. Neuroprotective and modulates inflammation.		
Taurodeoxycholic Acid	TDCA	Involved in bile acid metabolism and lipid metabolism. Modulates the gut microbiome and has anti-inflammatory properties.		
7-Ketolithocholic Acid	7-KLCA	Involved in bile acid metabolism and cholesterol metabolism. Modulates the gut microbiome.		
Dehydrolithocholic Acid	DHLCA	Involved in bile acid metabolism and modulation of the gut microbiome. Exerts anti-inflammatory properties.		
Hyodeoxycholic Acid	HDCA	May be negatively correlated with fecal calprotectin levels in UC. May have anti-atherosclerotic effects. Involved in lipid absorption and cholesterol metabolism. Modulates the gut microbiome.		
Alloisolithocholic Acid	Allolso-LCA	Positive associations with longevity and exerts antibacterial effects against gram-positive pathogens.		
3-Dehydrocholic Acid	3-DHCA	Involved in bile acid metabolism and exerts anti-inflammatory properties.		

^{*} On StoolOMX, the LCA value is a summation of LCA + Allo-LCA.



INDIVIDUAL SHORT CHAIN FATTY ACIDS RESULTS

Table 15.

SACCHAROLYTIC STRAIGHT CHAIN FATTY ACIDS (SCFA)			
MARKER	GUIDE		
Acetate	Most abundant SCFA. Involved in lipid synthesis and appetite regulation, maintains energy balance and metabolic homeostasis. Resists oxidation and mitochondrial stress.		
Butyrate	Primary energy source for colonic cells. Supports intestinal barrier integrity, reduces intestinal inflammation, promotes motility, enhances fatty acid oxidation, inhibits tumor cell progression, and fosters a balanced microbiome.		
Propionate	Supports intestinal barrier integrity, impacts energy balance, gluconeogenesis, and lipid metabolism. Involved in appetite regulation. Proposed as a biomarker for IBS.		
Valerate	Stimulates intestinal epithelium growth, inhibits colon cancer cell production, modulates glucose and lipid metabolism. Antimicrobial effects against <i>C. difficile</i> .		
Caproate	Antimicrobial effects against <i>C. difficile</i> .		

PROTEOLYTIC BRANCHED CHAIN FATTY ACIDS (BCFA)				
MARKER	GUIDE			
Isobutyrate	Stimulates colonic sodium absorption. Elevated levels in stool may be associated with cortisol levels and depression.			
Isovalerate	Elevated levels in stool may be associated with depression by its influence on the gut flora, metabolic pathways, and inflammatory pathways. It can also interfere with neurotransmitter release. There is an association between elevated levels of isovalerate and increased cortisol levels.			
2-Methylbutyrate	Branched chain fatty acid.			
Isocaproate	Produced primarily through the fermentation of branched chain amino acids (BCAAs), particularly leucine. Energy source for many tissues, especially during fasting or low carbohydrate intake.			

SCFAs are known to have numerous health effects and can enhance fecal excretion of bile acids.

Low levels of SCFAs are linked to various conditions, including irritable bowel syndrome, obesity, and inflammatory bowel disease.





Stool*OMX*™

Advanced Bile Acid Testing and Short-Chain Fatty Acid Evaluation in One GI-MAP Add-on Panel



GI-MAP PATTERNS

Understanding Common Dysbiosis Patterns With GI-MAP

INSUFFICIENCY DYSBIOSIS

Insufficiency dysbiosis is characterized by low levels of beneficial bacteria that provide critical support for healthy intestinal and immune function. Insufficient levels of beneficial bacteria may result in an elevated risk of intestinal infections, increased intestinal barrier permeability, decreased protective factors such as secretory IgA, and increased inflammation. Lack of keystone bacteria is common in autoimmune, allergic, and chronic inflammatory conditions.

Table 16.

Markers Characterizing Insufficiency Dysbiosis				
Commensal/Keystone Bacteria: (low levels)	Bacteroides fragilis Bifidobacterium spp. Enterococcus spp. Escherichia spp. Lactobacillus spp. Akkermansia muciniphilia Faecalbacterium prausnitzii Roseburia spp.			
Phyla Microbiota: (low levels)	Bacteroidetes Firmicutes			
Associated Intestinal Health Markers:	Secretory IgA (often low to very low levels) Zonulin (sometimes elevated)			





INFLAMMATORY DYSBIOSIS

Inflammatory dysbiosis is characterized by moderate to high levels of certain pathogens, normal microbiota, and opportunistic microbes that promote inflammation and increased intestinal permeability. Many pro-inflammatory microbes are gram-negative bacteria that belong to the Proteobacteria phylum and produce a form of lipopolysaccharide (*LPS*) that is a potent activator of inflammatory responses. This pattern is common in chronic immune and inflammatory conditions.

Table 17.

Markers Characterizing Inflammatory Dysbiosis				
Pathogens (low to high levels)	Campylobacter C. difficile Pathogenic E. coli Salmonella Vibrio cholerae Yersinia enterocolitica Giardia			
Commensal/Keystone Bacteria (high levels)	Escherichia spp. Enterobacter spp.			
Opportunistic Bacteria, Yeast, and Protozoa (moderate to high levels)	Morganella spp. Pseudomonas spp. Pseudomonas aeroginosa Desulfovibrio spp. Citrobacter spp. Citrobacter freundii Klebsiella spp. Klebsiella pneumoniae Proteus spp. Proteus mirabilis Fusobacterium spp. Candida spp. Candida albicans Parasitic protozoa (specifically Giardia and Blastocystis hominis)			
Associated Intestinal Health Markers:	β-Glucuronidase (may be elevated) Occult Blood-FIT (may be elevated) Secretory IgA (often low levels, but sometimes elevated) Calprotectin (often elevated, but sometimes very low levels) Eosinophil Activation Protein (EDN/EPX) (may be elevated) Zonulin (may be elevated in some cases)			



DIGESTIVE DYSFUNCTION DYSBIOSIS

Dysbiosis associated with digestive dysfunction is very common, and is often due to low stomach acid (hypochlorhydria), insufficient bile acids, poor digestion (pancreatic insufficiency or brush border enzyme deficiency), reduced absorption, and altered gastrointestinal motility. Altered digestion and motility can result in imbalances in the microbiome, characterized by overgrowth of certain species. Symptoms associated with digestive dysfunction include but are not limited to: excessive gas and bloating, abdominal discomfort, dyspepsia, heart burn, gastroesophageal reflux (GERD), constipation or diarrhea, food sensitivities and intolerances.

Table 18.

Markers Associate	ed with Digestive Dysfunction		
Pathogens (low to high levels)	Most types, especially if multiple pathogens are present		
H. pylori (moderate to high levels)	Helicobacter pylori (with or without virulence factors)		
Commensal/Keystone Bacteria (high levels)	Enterococcus Lactobacillus		
Phyla Microbiota (high levels)	Bacteroidetes and/or Firmicutes		
Opportunistic Bacteria, Yeast, and Protozoa (moderate to high levels)	Bacillus spp. Enterococcus faecalis Enterococcus faecium Morganella spp. Staphylococcus spp. Staphylococcus aureus Streptococcus spp. Methanobacteriaceae (family) Desulfovibrio spp. Klebsiella pneumoniae Prevotella Candida spp. Candida albicans Parasitic protozoa		
Intestinal Health Markers:	Elastase-1 (often low to moderately low levels) Steatocrit (sometimes elevated)		



PATTERNS ASSOCIATED WITH IMMUNE-MEDIATED FOOD REACTIONS

The GI-MAP® (Microbial Assay Plus) can be used in conjunction with IgG Food Explorer™ and IgE Allergy Explorer™ to dig deeper into the root cause of adverse food reactions. If the gut barrier is permeable and/or digestion is suboptimal, maldigested food proteins can trigger immune system responses.

These are the key patterns to look for on the GI-MAP that are connected food sensitivities and allergies.

Table 19.

Food Intolerand	ce, Allergy, and Adve	erse Food Re	actions Pattern
Histamine Producing	Morganella spp.		Opportunists (page 3)
	Pseudomonas spp.		
	Pseudomonas aeroginosa		
	Citrobacter freundii	High	
	Klebsiella spp.		
Bacteria	Klebsiella pneumoniae		
	Proteus spp.		
	Proteus mirabilis		
	Enterobacter spp.		
	Escherichia spp.		
	Fusobacterium spp.		
Mast Cell-Activating Bacteria	H. pylori	High	H. pylori Panel (page 2)
	Enterococcus faecalis	High	Opportunists (page 3)
	Pseudomonas aeruginosa		
	Staphylococcus aureus		
Bacteria	Streptococcus spp.		
	Candida spp.	High	Fungi/Yeast (page 3)
	Candida albicans		
	Lipopolysaccharide producers	High	Throughout report
Commensal/Keystone Bacteria	Lactobacillus spp.	Low	Commensal/Keystone Bacteria (page 2)
Intestinal Health Markers	Eosinophil Activation Protein (EDN/EPX)		Intestinal Health Markers (page 4)
	SIgA	High	
	Anti-gliadin IgA		
	Zonulin		

Continued...



Table 19a.

Gut Barrier Permeability ("Leaky Gut") Pattern				
Intestinal Permeability	Any Pathogen	High	Pathogens (page 1)	
	Lactobacillus spp.	Low	Commensal/Keystone	
	Akkermansia muciniphila	Low; <dl< td=""><td colspan="2">Bacteria (page 2)</td></dl<>	Bacteria (page 2)	
	Candida albicans	High	Fungi/Yeast (page 3)	
	Anti-gliadin IgA	High	Intestinal Health Markers (page 4)	
	Zonulin			
Low Butyrate/SCFA Production	Faecalibacterium prausnitzii	Low; <dl< td=""><td rowspan="3">Commensal/Keystone Bacteria (page 2)</td></dl<>	Commensal/Keystone Bacteria (page 2)	
	Roseburia spp.	Low		
	Firmicutes phylum			
Poor Mucosal Health	Bifidobacterium spp.	Low; <dl< td=""><td></td></dl<>		
	Escherichia spp.	1	Commensal/Keystone Bacteria (page 2)	
	Lactobacillus spp.	Low		
	Akkermansia muciniphila	Low; dl		
	Bacteroidetes phylum	Low		

Table 19b.

Digestive Insufficiency Pattern			
Digestive Insufficiency	Firmicutes phylum	High	Commensal/Keystone Bacteria (page 2)
	Bacteroidetes phylum		
	Enterococcus spp.		
	Lactobacillus spp.		
	Akkermansia muciniphila		
	Bacillus spp.	High	Opportunists (page 3)
	Enterococcus faecalis		
	Enterococcus faecium		
	Staphylococcus spp.		
	Staphylococcus aureus		
	Streptococcus spp.		
	Desulfovibrio spp.		
	Methanobacteriaceae (family)		
	Fusobacterium spp.		
	Steatocrit	Detected; High Intestinal	Intestinal Health Markers
	Pancreatic Elastase-1	Low	(page 4)





GAS & HISTAMINE PRODUCERS ON GI-MAP

Table 20.

Primary Hydrogen Producers

Faecalibacterium prausnitzii

Roseburia spp.

Bacteroidetes phyla

Firmicutes phyla

Primary Methane Producers Methanobacteriaceae (family)

Primary Hydrogen Sulfide Producers Escherichia spp.

Enterobacter spp.

Desulfovibrio spp.

Morganella spp.

Pseudomonas aeruginosa

Staphylococcus aureus

Citrobacter spp.

Citrobacter freundii

Klebsiella spp.

Klebsiella pneumoniae

Proteus spp.

Proteus mirabilis

Fusobacterium spp.

Bacteroides fragilis



MICROBE CATEGORIES

Table 21.

Primary Hydrogen Producers

Faecalibacterium prausnitzii

Roseburia spp.

Bacteroidetes phyla

Firmicutes phyla

Primary Methane Producers

Methanobacteriaceae (family)

Primary Hydrogen Sulfide Producers

Bacteroides fragilis

Escherichia spp.

Enterobacter spp.

Desulfovibrio spp.

Morganella spp.

Pseudomonas aeruginosa

Staphylococcus aureus

Citrobacter spp.

Citrobacter freundii

Klebsiella spp.

Klebsiella pneumoniae

Proteus spp.

Proteus mirabilis

Fusobacterium spp.

Histamine Producing Bacteria

Lactobacillus spp.

Morganella spp.

Pseudomonas

Pseudomonas aeruginosa

Citrobacter freundii

Klebsiella

Klebsiella pneumoniae

Proteus

Proteus mirabilis

Enterobacter spp.

Escherichia spp.

Fusobacterium spp.

Mast Cell-Activating Microbes

H. pylori

Enterococcus faecalis

Pseudomonas aeruginosa

Staphylococcus aureus

Streptococcus spp.

Candida spp.

Candida albicans

Lipopolysaccharide producers (see LPS list)

Lipopolysaccharide (LPS) Producing Bacteria

Escherichia spp.

Enterobacter spp.

Morganella spp.

Pseudomonas spp.

Pseudomonas aeruginosa

Citrobacter spp.

Citrobacter freundii

Klebsiella spp.

Klebsiella pneumoniae

Proteus

Proteus mirabilis





METHODOLOGY

Quantitative PCR (qPCR) or Real-Time Polymerase Chain Reaction (RT-PCR) — provides you with true quantitative values. It helps differentiate trace levels of an organism from frank elevations indicative of active infection.

SPECIMEN REQUIREMENTS

Single Stool Sample — at ambient room temperature in specimen vial provided.

TEST ORDERING OPTIONS

- GI-MAP® | GI Microbial Assay Plus
- GI-MAP® + Zonunlin
- · GI Pathogens Profile
- H. pylori Profile
- Zonulin Profile (Stool)
- Calprotectin (Stand-alone)
- Gluten Peptide (Add-on or stand-alone)
- Universal Antibiotic Resistance Genes Panel (Add-on only)
- StoolOMX[™] (Add-on only)

Build Your GI-MAP Knowledge With Our Complimentary Educational Services at DSL Academy:

DSL Academy (DSLA) is an exclusive, online learning platform designed to help you better interpret the GI-MAP test. The DSLA curriculum is taught by the DSL Medical Education Team and is available to all active DSL account holders. *Learn more at: www.diagnosticsolutionslab.com/clinicians/dsl-academy.*





Diagnostic Solutions laboratory

Phone: 877-485-5336

Email: cs@diagnosticsolutionslab.com Web: diagnosticsolutionslab.com


SAMPLE REPORT



Scan with Your Camera

Download a GI-MAP® Sample Report

diagnosticsolutionslab.com/tests/gi-map



INTERPRETIVE GUIDE

GI-MAP® – Unparalleled DNA Based Stool Testing

Our mission: to deliver innovative, accurate and clinically relevant diagnostic testing in a timely and cost-effective manner



GI Microbial Assay Plus

